
Mouse Anti Human CD19 FITC/CD5 PE

Catalog # CD19-D3-100 **Size** 100 tests
Catalog # CD19-D3-25 **Size** 25 tests

Clone Number: HIB19/HISM2
Isotype: IgG1/IgG1

Product Form: Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.