

<b>Catalog#</b>	<b>CD45-D1-100</b>	<b>Size</b>	<b>100 tests</b>
<b>Catalog#</b>	<b>CD45-D1-25</b>	<b>Size</b>	<b>25 tests</b>

---

**Mouse Anti Human CD45 FITC/CD14 PE**

**Clone Number:** HI30/MEM-15  
**Isotype:** IgG1/IgG1  
**Volume/Quantity:** 25T

**Product Form:** Dual color combination consisting of FITC conjugated and R-PE conjugated monoclonal antibodies mixed in optimal ratio.

**USAGE**

The reagent is tested for flow cytometric analysis using 20µl/10<sup>6</sup> cells or 100µl peripheral blood cells.

**STORAGE**

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.