

Catalog# **CD45-T4-100** **Size** **50 tests**
Catalog# **CD45-T4-25** **Size** **25 tests**

Mouse Anti Human CD45 FITC/CD8 PE/CD3 PE-Cy5

Clone Number: HI30/HIT8a/HIT3b
Isotype: IgG1/IgG1/IgG1
Volume/Quantity: 50T

Product Form: Triple color combination consisting of FITC, R-PE and PE-Cy5 conjugated monoclonal antibodies mixed in optimal ratio.

USAGE

The reagent is tested for flow cytometric analysis using 20µl/10⁶ cells or 100µl peripheral blood cells.

STORAGE

Storage at 4°C, should not be frozen and avoid prolonged exposure to light.

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400
Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com