

<b>Catalog#</b>	<b>CD58P-100</b>	<b>Size</b>	<b>100 tests</b>
<b>Catalog#</b>	<b>CD58P-25</b>	<b>Size</b>	<b>25 tests</b>

---

**Mouse Anti-Human CD58-PE conjugate**

**PRODUCT INFORMATION**

**CLONE:** HI58a  
**ISOTYPE:** Mouse IgG1  
**Product Forms:** Purified, FITC conjugation, PE conjugation.  
**SIZE:** 25 Tests

**DESCRIPTION**

CD58 McAb recognizes a 55-70KD membrane glycoprotein known as the lymphocyte function associated antigen-3 (LFA-3). CD58 antigen has wide distribution over many tissues, both hematopoietic and non-hematopoietic cells, including leucocytes, erythrocytes, endothelial cells, epithelial cells and fibroblasts. There are two isoforms of CD58 antigen: a glycosylphosphatidylinositol (GPI)- anchor form and a transmembrane form. Both forms are expressed on the same cell type, and however, erythrocytes only have the GPI- anchored isoform. CD58 is the ligand for CD2, and functions as an adhesion molecule and mediates its signal through CD2, which is expressed on T cells and NK cells.

**PREPARATION**

The monoclonal antibody is purified from ascites by protein G affinity chromatography.

**USAGE**

The purified reagent is effective for indirect immunofluorescence staining of human cells for flow cytometric analysis.

The conjugated reagent (FITC, R-PE) is tested for flow cytometric analysis using 20µl/10<sup>6</sup> cells or 100µl peripheral blood cells.

**STORAGE**

For purified forms, long term storage at -20°C.

For conjugated forms, storage at 4°C, should not be frozen and avoid prolonged exposure to light.

---

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400  
Fax: +91-11-42208444 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)