



India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208444 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com

Item Name : Chicken Adult Tissues Protein lysate, Universal

Cat# CHUL-4567

Size: 100 ug

Species: Chicken

Tissue type Chicken Normal Tissue, protein pools from multiple donors

Donor Information: Male: 8 month(s) old. Pool of 1 donor(s)
Female: 8 month(s) old. Pool of 1 donor(s)

Buffer:

Total tissue protein lysate was prepared in Hepes buffer, pH 7.9, KCl, dexoycholate, NP-40%, EDTA, sucrose, and a cocktail of protease inhibitors. Lysate is typically supplied at a concentration of 2 mg/ml (or see lot specific conc. on the vial).

Storage and stability:

Store at 2-8°C for short term use. Long term storage at -70°C in suitable aliquots. Avoid repeated freeze thaw cycles. Under proper storage conditions the shelf life ~6 months from the date of receipt.

Recommended Usage

Western Blot (10 µg to 20 µg per lane is recommended for a mini gel). Immunoprecipitation, Electrophoresis, Enzymatic activity analysis, Protein-protein interaction, tissue specific expression.

All items are for in vitro research use only.

Related items

#90100	Western blot recycling kit
CHUL-4567	Chicken Adult Tissues Protein lysate, Universal
CHUC-565	Chicken (Normal/Adult) Universal (pool of tissue cDNA)
DOUL-4568	Dog Adult Tissues Protein lysate, Universal
DOUC-566	Dog (Normal/Adult) Universal (pool of tissue cDNA)
ERA12-C	Human Cell-18 protein lysate for WB positive control for ERA12-A protein (PR) IgG #2
HSP703-C	Human HeLa cells (normal/non-heat shocked) total protein lysate for western control for hsp
HSP704-C	Human HeLa cells (Heat shocked) total protein lysate for western control for hsp
HUUL-4565	Human Adult Tissues Protein lysate, Universal
MKUL-4569	Monkey Adult Tissues Protein lysate, Universal
MOUL-4566	Mouse Adult Tissues Protein lysate, Universal
PR12-C	Human Cell-30 protein lysate for WB positive control for PR12-M r protein (PR) IgG #2
RTUL-4570	Rat Adult Tissues Protein lysate, Universal

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

CHUL-4567 150313P

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com