

Product Specification Sheet

**Chloride Channel-1 (CLC-1 or CLCN1) Antibodies**

Cat. # CLC12-P	Rat CLC1 Control Peptide #2	<b>SIZE:</b> 100 ug
Cat. # CLC12-S	Rabbit Anti-rat CLC-1 antiserum #2	<b>SIZE:</b> 100 ul
Cat. # CLC12-A	Rabbit Anti-rat CLC-1 Ig G #2 (affinity pure)	<b>SIZE:</b> 100 ug

Chloride is a critical component of all living cells. Voltage-gated chloride channels regulate cellular traffic of chloride ion. The chloride channels (CIC or CLC) performs several functions including the regulation of cell volume, membrane potential stabilization, signal transduction, and transepithelial transport. Mutations in CIC genes have been linked with several human diseases including myotonias (Thomsen's disease), cystic fibrosis, Bartters syndrome type III, Dent's disease, and X-linked recessive nephrolithiasis. In mammals, CLC proteins form a superfamily of at least 9 different genes (CLC1-7 also known as CLCN1-7 and CLK1-2 or CLCKa and CLCKb). Additional forms of these proteins are obtained by alternative splicing. All CLC proteins (~700-1000 aa) are predicted to contain 10 (possibly 12) transmembrane domains. Except CLC-1 and CLC-K1/K2 that are specific for kidney, most other CLC are widely distributed in various tissues.

Rat CLC-1 is 994 aa membrane protein (human CLC-1 988 aa) (1). It is predominantly expressed in skeletal muscles. Defects in CLC1 (CLCN1) are the cause of autosomal recessive generalized myotonia (Becker's disease) (RGM) and autosomal dominant myotonia congenita (Thomsen's disease; MC) which are characterized by skeletal muscle stiffness (delayed relaxation) that is a result of muscle membrane hyperexcitability

**Source of Antigen and Antibodies**

<b>Antigen</b>	19aa peptide of rat CLC1; <b>Designated (CLC12-P or control peptide)</b> . epitope location ~ N-terminus
<b>Ab Host/type</b>	Rabbit, polyclonal Unpurified antiserum (cat #CLC12-S) Aff pure IgG1 (cat #CLC12-A) purified over antigen-agarose column
<b>2-ab</b>	<b>Goat Anti-rabbit IgG-HRP</b> cat # 20320 (AP, biotin, FITC conjugates also available)
<b>-ve control IgG</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied 0.05% azide, **Reconstitute** powder in 100 ul PBS

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder** in PBS at 1mg/ml

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at -20OC and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure antibody using ECL technique). See refs in 2

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

**Histochemistry & Immunofluorescence:** We recommend the use of affinity purified antibody at 1-20 ug/ml in paraformaldehyde fixed sections of tissues. See refs in 2

**Specificity & Cross-reactivity**

The 19 AA rat CLC12-P control peptide 94% conserved in mouse, 84% in human and 57% in dog CLC-1. No significant sequence homology is detected with other CLCs or other proteins. Antibody cross-reactivity in various species has not been studied. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at the web site).

**General References:** Steinmeyer K et al (1991) Nature 354, 301; Koch MC et al (1992) Science 257, 797; George AL et al (1993) Nature Genet. 3, 305; Steinmeyer K et al (1994) EMBO J 13, 737; Lorenz C et al (1994) Human Mol. Genet. 3, 941; Heine R et al (1994) Hum. Mol. Genet. 3, 1123; George Al et al (1994) Hum. Mol. Genet. 3, 2071; Meyer-Kleine C et al (1995) Am. J. Hum. Genet. 57, 1325; Lehman Horn F et al (1995) Hum. Mol. Genet 4, 1397; Mailander V (1996) Am. J. Hum. Genet. 58, 317

**2. Citations for ADI Antibodies** (see updates at the web site)

Wheeler TM, 2007, J. Clin. Invest., 117: 3952 - 3957. IF  
Kanadia RN, 2006, PNAS, 103: 11748 - 11753 WB IHC  
Yu Z, 2006, J. Clin. Invest., Oct 2006; 116: 2663 – 2672, IF  
\*This product is for In vitro research use only.

**Related material available from ADI**

Antibodies CLC1-7 and CLC-K1; KCCL1-3; AQP-9 and RUT; OCT and OAT, AE-3, and NACX

CLC12-S-A-P 71217A

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400

Fax: +91-11-42208444 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)