

Product Specification Sheet

**Carboxyl Terminal Modulator Protein (CTMP) Antibodies**

<b>Cat.</b> CTMP11-S	Rabbit Anti-Human CTMP Antiserum # 1	<b>SIZE:</b> 100 ul
<b>Cat.</b> CTMP11-A	Rabbit Anti-Human CTMP IgG # 1(Aff Pure)	<b>SIZE:</b> 100 ug
<b>Cat.</b> CTMP11-P	Human CTMP Control peptide # 1	<b>SIZE:</b> 100 ug
<b>Cat.</b> CTMP11-C	Recombinant purified Human CTMP protein Control for WB	<b>SIZE:</b> 100 ul

Protein kinase B (PKB) is a major downstream target of receptor tyrosine kinases that signal via phosphatidylinositol 3-kinase (PIK3). Upon cell stimulation, PKB is translocated to the plasma membrane, where it is **phosphorylated on thr308** in the catalytic domain and **ser473** in the C-terminal regulatory domain. Several protein partners for PKB have been identified. **CTMP, a carboxyl terminal modulator protein**, binds to PKB- $\alpha$  and negatively regulates PKB activity by inhibiting phosphorylation of these PKB residues. CTMP (mouse 248-aa, human 240-aa, chromosome 1q21, **22-26 kDa**) predominantly expressed in skeletal muscle, testis, uterus, brain, and kidney, with lower levels in heart, liver, and lung. The presence of multiple transcripts suggested that CTMP undergoes alternative splicing.

**Sources of antigen and antibodies**

<b>Antigen</b>	17aa peptide of Human <b>CTMP (1); Designated (CTMP11-P or Control peptide /blocking peptide)</b> conjugated to KLH. Epitope location ~ C-terminus
<b>Ab Host/type</b>	Rabbit, Polyclonal unpurified antiserum ( <b>#CTMP11-S</b> ) and IgG, purified over antigen-agarose (Cat # <b>CTMP11-A</b> )
<b>2-Ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

Human recombinant CTMP protein (protein accession #Q96KR2, 240aa ~27 kda ) was expressed as GST-tag protein, GST-tag removed, and purified (>95%). For Western blot +ve control (**Cat # CTMP11-C**) is supplied in SDS-PAGE sample buffer (reduced). Load 10 ul/lane of **CTMP11-C** for good visibility with antibody Cat # **CTMP11-S**. Store at -20oC in suitable size aliquots. SDS may crystallize in cold conditions. It should redissolve by warming before taking it from the stock. It should be heated once prior to loading on gels. If the product has been stored for several weeks, then it may be preferable to add 5 ul of fresh 2x sample buffer per 10 ul of the **CTMP11-C** solution prior to heating and loading on gels. This preparation is not biologically active. It is not suitable for ELISA or other applications where native protein is required. Do not freeze, thaw, or heat repeatedly

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**  
100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute** powder in 100 ul PBS

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in Buffer: PBS+0.1% BSA  
**Reconstitute powder** in PBS at 1mg/ml

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at 20°C and powder at 4°C or -20°C..

**Long-term:** at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20°C or below.

**Shipping:** 4oC for solutions and room temp for powder.

**Recommended Usage**

**Western Blotting** (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure using Chemiluminescence technique).

**ELISA** (1:10K-1:100K; using 50-100 ng of control peptide/well).

**Histochemistry & Immunofluorescence:** not tested. We recommend the use of affinity pure antibody at 2-20 ug/ml.

**Specificity & Cross-reactivity**

The Human CTMP11-P control peptide sequence is 100% conserved in mouse CTMP. Antibody crossreactivity in various species is not known. The control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at: the web site).

**General References:** (1) Maria S-M et al (2001) Science 294, 374-380.

\*This product is for in vitro research use only.

**Related material available from ADI:**

Antibodies against AKT1, AKT2, AKT3

**Western Blot recycling kit** (Use the same blot to probe with multiple antibodies AKT1, AKT or AKT3,) **recycle blot at room temp in 5-10 min;** No mercaptoethanol or heating required).

CTMP11-S-A-P-C

71210S

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A;

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)