

Product Specification Sheet

**Connexin Mimetic Peptides**

Cat. # Cx2409-P-1	Human/Mouse/Rat <sup>32</sup> Cx-Gap <sub>24</sub> Control Peptide	<b>SIZE:</b> 1 mg
Cat. # Cx2409-P-5	Human/Mouse/Rat <sup>32</sup> Cx-Gap <sub>24</sub> Control Peptide	<b>SIZE:</b> 5 mg
Cat. # Cx2410-PS-1	Human/Mouse/Rat <sup>32</sup> Cx-Sc-Gap <sub>24</sub> Scrambled Peptide	<b>SIZE:</b> 1 mg
Cat. # Cx2410-PS-5	Human/Mouse/Rat <sup>32</sup> Cx-Sc-Gap <sub>24</sub> Scrambled Peptide	<b>SIZE:</b> 5 mg

**Gap junctions** contain channels that allow the passage of ions and small molecules between adjacent cells. This intercellular communication has been implicated in the coordination of cellular responses to intracellular signaling molecules.

**Connexins** are a family of integral membrane proteins that oligomerize to form intercellular channels that are clustered at gap junctions. Connexins share a common transmembrane (TM) topology. Each connexin is thought to contain 4 TM domains, with two extracellular and three cytoplasmic regions.

**Connexin-mimetic peptides**, i.e. peptides that are identical with a short amino acid sequence on the Cx subunit have recently emerged as tools to block Connexin hemi-channels (CxHcs) with little or no immediate effects on gap junctions.

**Gap 24** corresponds to amino acid sequence on the intracellular loop of Cx32. Gap24 inhibits transfer of fluorescent dyes, electrical coupling and synchronized Ca<sup>2+</sup> oscillations in smooth-muscle cells by inhibiting CXhc32 and not gap junction (1).

**Source of Peptide**

<b>Cat #s</b>	<ul style="list-style-type: none"> <li>Cx2409-P-1</li> <li>Cx2409-P-5</li> </ul>
<b>Sequence</b>	G-H-G-D-P-L-H-L-E-E-V-K
<b>Identity</b>	Human/Mouse/Rat <sup>32</sup> Cx-Gap <sub>24</sub> <b>Control Peptide</b>
<b>Location</b>	110 - 121 aa
<b>Purity</b>	>90%; HPLC Purified
<b>Mol. Wt.</b>	1330.46
<b>Form</b>	Lyophilized Powder
<b>Solubility</b>	Soluble in DMF

<b>Cat #</b>	<ul style="list-style-type: none"> <li>Cx2410-PS-1</li> <li>Cx2410-PS-5</li> </ul>
<b>Sequence</b>	L-P-D-G-V-K-H-G-E-E-L-H
<b>Identity</b>	Human/Mouse/Rat <sup>32</sup> Cx-Gap <sub>24</sub> <b>Scrambled Peptide</b>
<b>Location</b>	110 - 121 aa
<b>Purity</b>	>90%; HPLC Purified
<b>Mol. Wt.</b>	1330.46
<b>Form</b>	Lyophilized Powder
<b>Solubility</b>	Soluble in DMF

**Storage**

**Short-term:** unopened, undiluted liquid vials for less than a week at 4°C.

**Long-term:** Store at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw or store working, diluted solutions.

**Stability:** 6-12 months at -20°C or below.

**Shipping:** 4°C for solutions and room temp for powder.

**Recommended Usage**

**General References:**

1) De Vuyst E, et al., (2006) EMBO J. 25: 34-44.

\*This product is for In vitro research use only.

**Related material available from ADI**

Control and scrambled Peptides to mouse, rat and human Gap26, Gap27 and CxHc-32

Cx2409-10-P-1-5

71218J