

Petri Dish, 100 x 20 mm, Slippable

Product ID D2483

Description

Packed in cases of 400 each with 20 sleeves of 20 dishes. Sterile.

- Sterile, 100 x 20 mm Petri dish
- Slippable
- ISO Mark target for automation
- Disposable, one-time use (not autoclavable)

Product Information

Autoclavability	Not autoclavable
Diameter	100 mm
Material	Polystyrene
Product Number	D2483
Height	20 mm

Shipping and Storage

Tariff Code	3923.10.0000
-------------	--------------

PhytoTechnology Laboratories®

14610 W 106th St. Lenexa, KS 66215

India Contact:

Life Technologies (India) Pvt. Ltd.

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com