

Enzyme Conjugate for Potato virus Y-common strain (PVY-c)



Test Label: Alkaline Phosphatase

Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

Option

Enzyme conjugate for Potato virus Y-common strain (PVY-c), alkphos, 500 testwells
ECA 20700/0500

Enzyme conjugate for Potato virus Y-common strain (PVY-c), alkphos, 1000 testwells
ECA 20700/1000

Enzyme conjugate for Potato virus Y-common strain (PVY-c), alkphos, 5000 testwells
ECA 20700/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com