

## Enzyme Conjugate for Beet western yellows virus (BWYV)



**Test Label:** Alkaline Phosphatase

**Test Format:** Triple Antibody Sandwich (TAS) ELISA

### Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

### Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

### Option

Enzyme conjugate for Beet western yellows virus (BWYV), alkphos, 500 testwells  
ECA 32100/0500

Enzyme conjugate for Beet western yellows virus (BWYV), alkphos, 1000 testwells  
ECA 32100/1000

Enzyme conjugate for Beet western yellows virus (BWYV), alkphos, 5000 testwells  
ECA 32100/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)