

## Enzyme Conjugate for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt)



**Test Label:** Alkaline Phosphatase  
**Test Format:** Double Antibody Sandwich (DAS) ELISA  
**Intended Use:**

This enzyme conjugate is intended for use with its corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

### Includes:

- Alkaline phosphatase enzyme conjugate
- [User Guide](#)

### Option

Enzyme conjugate for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), alkphos, 500 testwells  
ECA 44100/0500

Enzyme conjugate for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), alkphos, 1000 testwells  
ECA 44100/1000

Enzyme conjugate for *Clavibacter michiganensis* subsp. *tessellarius* (Cmt), alkphos, 5000 testwells  
ECA 44100/5000

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)