

Enzyme Conjugate for Watermelon silver mottle virus & Groundnut bud necrosis virus (WSMoV & GBNV)



Test Label: Alkaline Phosphatase
Test Format: Triple Antibody Sandwich (TAS) ELISA

Intended Use:

This detection antibody (bottle A) and enzyme conjugate (bottle B) are intended for use with their corresponding capture antibody (CAB).

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Detection antibody, bottle A
- Alkaline phosphatase enzyme conjugate, bottle B
- [User Guide](#)

Option

Enzyme conjugate for Watermelon silver mottle virus & Groundnut bud necrosis virus (WSMoV & GBNV), alkphos, 500 testwells
ECA 61500/0500

Enzyme conjugate for Watermelon silver mottle virus & Groundnut bud necrosis virus (WSMoV & GBNV), alkphos, 1000 testwells
ECA 61500/1000

Enzyme conjugate for Watermelon silver mottle virus & Groundnut bud necrosis virus (WSMoV & GBNV), alkphos, 5000 testwells
ECA 61500/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com