

Anti-animal Enzyme Conjugate for *Xanthomonas campestris* pv. *armoraciae* (Xcarm)



Test Label: Alkaline Phosphatase
Test Format: Antigen coated plate (ACP) ELISA

Intended Use:

This anti-animal enzyme conjugate is intended for use with its corresponding detection antibody in ACP ELISA format.

Test performance using other capture antiserum has not been evaluated and is not guaranteed.

Includes:

- Alkaline phosphatase enzyme conjugate
- [User Guide](#)

Option
Xcarm enzyme conjugate, alkphos, 500 testwells ECA 88800/0500
Xcarm enzyme conjugate, alkphos, 1000 testwells ECA 88800/1000
Xcarm enzyme conjugate, alkphos, 5000 testwells ECA 88800/5000

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com