


Product Specification Sheet

Recombinant Zaire Ebola virus secreted glycoprotein (sGP)

- **Cat # EVGP23-R-10** Recombinant (HEK) Zaire Ebola virus secreted glycoprotein (sGP, Mayinga/76, 33--364aa, Flag-myc-his-tag, >95%), low endotoxin **SIZE:** 10 ug
- **Cat # EVGP23-R-100** Recombinant (HEK) Zaire Ebola virus secreted glycoprotein (sGP, Mayinga/76, 33--364aa, Flag-myc-his-tag, >95%), low endotoxin **SIZE:** 100 ug

Ebola virus (EBOV, formerly Zaire ebolavirus) causes severe disease in humans and in nonhuman primates in the form of viral hemorrhagic fever. Zaire ebolavirus is a virological taxon included in the genus Ebolavirus, family Filoviridae, order Mononegavirales. The species has a single virus member, Ebola virus (EBOV). **Ebolavirus species Zaire (ZEBOV)** causes highly lethal hemorrhagic fever, resulting in the death of **90%** of patients within days. Most information on immune responses to ZEBOV comes from in vitro studies and animal models. Ebola Zaire attacks every organ and tissue in the human body except skeletal muscle and bone. Ebola is classified as a **Level 4** pathogen (higher than AIDS) with a 2 to 21 days (7 to 14 days average) incubation period. There are currently four known strains of Ebola: **Zaire, Sudan, Reston and Tai**. All cause illness in sub-human primates. Only Ebola Reston does not cause illness in humans. The mortality rate of Ebola victims is between 60% and 90%; with Ebola Sudan at 60% and Ebola Zaire at 90%.

The virions are tubular in general form but variable in overall shape and may appear as the classic shepherd's crook or eyebolt. Ebola virions consist of 7 structural proteins. At the center is the **helical ribonucleocapsid**, which consists of the genomic RNA wrapped around a polymer of **nucleoproteins (NP)**. Associated with the ribonucleoprotein is the RNA-dependent **RNA polymerase (L)** with the **polymerase cofactor (VP35)** and a **transcription activator (VP30)**. The ribonucleoprotein is embedded in a matrix, formed by the major (VP40) and minor (VP24) matrix proteins. They are surrounded by a **lipid membrane** derived from the host cell membrane. The membrane anchors a glycoprotein (GP1,2) that projects 7 to 10 nm spikes away from its surface. While nearly identical to **Marburg virions** in structure, ebola virions are antigenically distinct.



The most common diagnostic methods are RT-PCR in conjunction with antigen-capture ELISA which can be performed in field or mobile hospitals and laboratories. There are currently no FDA-approved vaccines for the prevention of EVD. The most promising ones are DNA vaccines or are based on adenoviruses, vesicular **stomatitis Indiana virus (VSIV) or filovirus-like particles (VLPs)** as all of these candidates could protect nonhuman primates from Ebola virus-induced disease. DNA vaccines, adenovirus-based vaccines, and VSIV-based vaccines have entered clinical trials.

Source of Antigen

Recombinant Zaire Ebola virus sGP/pre-small/secreted glycoprotein (NP_066246, 331le-3641le) was expressed in HEK cells as a Flag-Myc-his-tag fusion protein (>95%). Purified protein is supplied in 30 mM HEPES, 150 mM NaCl, pH 7.0 (see lot sp. Conc. on the vial) or lyophilized in the same buffer. Store powder at -20oC. Reconstitute powder in water and store at -20oC in suitable size aliquots.

It is suitable for ELISA, Western or other applications where native protein is required. Do not freeze, thaw, or heat repeatedly.

Endotoxin: measures at <0.1 EU/ µg protein as determined by the LAL method.

Stability: 6-12 months at -20oC or below.

Shipping: 4oC for solutions and room temp for powder.

Recommended Usage

Western Blotting: load 20-100 ng/well.

ELISA (50-100 ng antigen/well).

Specificity

Zaire pre-small/sGP/small non-structural GP (Mayinga/76/33-364aa) is 100 % conserved in Ebola Gabon (1994), Zaire Ebola (1995), 75% in Tai-Forest sGP (1994) & Bundibugyo sGP, 66% in Reston (1989) sGPs.

References: Thomas W (2010) Archives of Virology 155 (12): 2083–103. Taylor D (2010) BMC Evolutionary Biology 10: 193. Feldmann H (2005) . A. Virus Taxonomy—Eighth Report of the International Committee on Taxonomy of Viruses. 645–653.

*This product is for In vitro research use only.

http://www.4adi.com/objects/catalog/product/extras/Ebola_Marburg_Vaccines_ELISA_Flr.pdf

- EVGP15-A Anti-Zaire Ebola virus glycoprotein (GP, 1-676aa/DNA vaccine) IgG,
- EVGP16-A Anti-Zaire Ebola virus glycoprotein (GP 1-652aa/DNA vaccine) IgG,
- EVGP16-R-10 Rec. (sf9) Sudan-Ebola virus glycoprotein (minus transmembrane domain, his-tag, 68 kda), purified
- EVGP17-R-10 Recombinant (sf9) Zaire-Ebola virus glycoprotein (minus transmembrane domain, his-tag, 68 kda), purified
- EVGP20-R-10 Recombinant (sf9) Zaire Ebola virus glycoprotein (GIN/2014/Kissidougou-C15, 1-650aa, his-tag at CT, >95%), low endotoxin
- EVGP21-R-10 Recombinant (HEK) Zaire Ebola virus glycoprotein (GIN/2014/Kissidougou-C15, 1-650aa, his-tag at CT, >95%), Low endotoxin
- BVRB11-R-10 Recombinant (HEK) Bundibugyo Ebola virus glycoprotein RBD domain (Uganda 2007, 54-201aa, Fc-tag at CT, >95%, low endotoxin)
- EVRB14-R-10 Recombinant (HEK) Zaire Ebola virus glycoprotein RBD domain (Mayinga 1976, 1-308 aa, his tag, >95%, low endotoxin)
- SVRB13-R-10 Recombinant (HEK) Sudan-Ebola virus RBD domain (Gulu, 1-320aa, Fc-tag at CT, >95%, low endotoxin)
- EVRB13-R-10 Recombinant (HEK) Zaire Ebola virus glycoprotein RBD domain (Mayinga 1976, 1-308 aa, Fc tag at CT, >95%, low endotoxin) purified
- SVGP24-R-10 Recombinant (HEK) Sudan-Ebola virus glycoprotein (Gulu, 1-637aa, his-tag at CT, >95% low endotoxin)
- EVNP13-A Anti-Zaire-Ebola virus nucleoprotein (EBOV NP, 1-739/DNA vaccine) IgG,
- EVP401-A Anti-Zaire-Ebola virus VP40 peptide (EBOV VP40) IgG,
- EVP401-C Rec. Zaire-Ebola virus VP40 protein control for Western
- EVP405-R-10 Rec. (E.coli) Zaire-Ebola virus VP40
- AE-320520-1 Human Anti-Ebola virus Nucleoprotein (NP) IgG ELISA Kit
- AE-320530-1 Human Anti-Ebola virus Nucleoprotein (NP) IgM ELISA Kit
- AE-320620-1 Human Anti-Zaire-Ebola virus glycoprotein (GP) IgG ELISA Kit
- AE-320720-1 Human Anti-Zaire-Ebola virus VP40 IgG ELISA Kit,
- AE-320730-1 Human Anti-Zaire-Ebola virus VP40 IgM ELISA Kit,
- AE-320800-48 Human Zaire-Ebola Virus antigen ELISA Kit, 48 tests,
- AE-320800-96 Human Zaire-Ebola Virus antigen ELISA Kit, 96 tests,

EVGP23-R-Recombinant-Zaire-Ebola-sGP-Protein 150421A

Alpha Diagnostic Intl Inc., 6203 Woodlake Center Dr, S an Antonio, T X 7 8 24 4 , U S A ;

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400 Fax: +91-11-42208444 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com