

# Glove, Autoclave, 18 in

Product ID G089

## Description

Beige cotton terry cloth autoclave gloves. For use up to 232 °C. 18 inch length.

## Product Information

Material	Cotton Terry Cloth
Product Number	G089
Length	18 in (45.72 cm)

## Shipping and Storage

Storage Temp.	Room Temperature
Tariff Code	6216.00.3800

**PhytoTechnology Laboratories®**

14610 W 106th St. Lenexa, KS 66215

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)