

Product Specification Sheet

**Glutamate Aspartate Transporter (GLAST/EAAT1) Antibodies**

<b>Cat #</b> GLAST11-P	Rat GLAST control peptide	<b>SIZE:</b> 100 ug
<b>Cat #</b> GLAST11-S	Rabbit Anti-Rat GLAST antiserum	<b>SIZE:</b> 100 ul
<b>Cat #</b> GLAST11-A	Rabbit Anti-Rat GLAST IgG, aff pure	<b>SIZE:</b> 100 ug

Glutamate is the main excitatory neurotransmitter in the brain. To date five glutamate Transporters have been cloned: **GLAST (EAAT1)**, **GLT1 (EAAT2)**, **EAAC1 (EAAT3)**, **EAAT4**, and **EAAT5**. These transporters are believed to be critical in reducing potentially toxic extracellular concentration of glutamate by rapid uptake into nerve terminals and glial cells. Glutamate transporters (525-573 AA) display about 55% homology and are predicted to contain up to 6-10 transmembrane domains. Immunolocalization studies indicate that **GLT1** is localized in astroglial cells throughout the brain and spinal cord. **EAAC1** is specific for certain neurons and purkinje cells, and specifically enriched in cortex, hippocampus, and caudate-putamen and confined to presynaptic and postsynaptic elements. It is also expressed in kidney, heart, lung and muscle. **GLAST** has been observed in both neuron and astroglia. It is most abundant in Bergmann glia, cortex, hippocampus and cerebellum. **EAAT4** has properties of ligand gated Cl-channel. It is localized mainly in cerebellar Purkinje cells in rat and human CNS. **EAAT5** has only been cloned from human. It is primarily expressed in retina.

**Source of Antigen and Antibodies**

<b>Antigen</b>	15-aa peptide from <b>rat GLAST (1)</b> ; <b>Designation (#GLAST11-P, control/blocking peptide)</b> conjugated to KLH. <b>Epitope location</b> ~ C-terminus
<b>Ab Host/type</b>	Rabbit, Polyclonal unpurified antiserum ( <b>#GLAST11-S</b> ) and IgG, purified over antigen-agarose (Cat # <b>GLAST11-A</b> )
<b>2-Ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as –ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute powder in 100 ul PBS**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder in PBS at 1mg/ml**

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at 20°C and powder at 4°C or -20°C..

**Long-term:** at –20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at –20°C or below.

**Shipping:** 4°C for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** 1:1K-5K for antiserum and 1-10 ug/ml for affinity pure IgG using Chemiluminescence technique. This antibody has been used in many studies (2).

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

**Histochemistry & Immunofluorescence:** We recommend the use of affinity purified IgG at 2-10 ug/ml in paraformaldehyde fixed sections of tissues (see published refs Deitch JS et al (2002).

**Specificity & Cross-reactivity**

The rat GLAST11-P peptide is 100% conserved in mouse, human, bovine, and 93% in chicken, dog. The GLAST11-P has no significant sequence homology to other glutamate transporters. Antibody crossreactivity in various species is not established. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at:the web site).

**General References: (1).** Tanaka, K. (1993) Neurosci. Lett. **159**:183; Storck, T. et al (1992) PNAS **89**:10955; Tanaka, K. (1993) Neurosci. Lett. **16**:149; Shashidharan, P. et al (1993) BBA **1216**:161; Rothstein, J. D. et al (1994) Neuron **13**:713

**(2) Citations of ADI's Antibodies** (see web site for updated list)

**(2).** Alexander GM et al (2000) J. Neurochem. 74, 1666; Ueda Y et al (2000) Exp. Brain Res. 133, 334-340; Ueda Y et al (2000) Epilepsy Res. 39, 201; Ueda Y et al (2001) J. Neurochem. 76, 892; Deitch JS et al (2002) J. Neurosci. 193, 117;

*\*This product is for In vitro research use only.*

**Related material available from ADI**

Anti-GLAST, EAAC1, GLT1, EAAT4, EAAT5, GTRAP41, GTRAP48, VGLUT1/BNPI, VGLUT-2/DNPI & GABA Transporters (GAT1-3)

GLAST11-S-A-P

71226S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)