

Product Specification Sheet

**Glutamate Transporter (GLT1/EAAT2) Antibodies**

<b>Cat #</b> GLT11-P	Rat GLT1 control peptide	<b>SIZE:</b> 100 ug
<b>Cat #</b> GLT11-S	Rabbit Anti-Rat Rat GLT1 antiserum	<b>SIZE:</b> 100 ul
<b>Cat #</b> GLT11-A	Rabbit Anti-Rat Rat GLT1 IgG, aff pure	<b>SIZE:</b> 100 ug

Glutamate is the main excitatory neurotransmitter in the brain. To date five glutamate Transporters have been cloned: **GLAST (EAAT1), GLT1 (EAAT2), EAAC1 (EAAT3), EAAT4, and EAAT5**. These transporters are believed to be critical in reducing potentially toxic extracellular concentration of glutamate by rapid uptake into nerve terminals and glial cells. Glutamate transporters (525-573 AA) display about 55% homology and are predicted to contain up to 6-10 transmembrane domains. GLT1 (also known as Excitatory Amino acid Transporter 2, EAT2; or Sodium-dependent glutamate/Aspartate transporter 2) is a glutamate transporter. It is a 573 amino acid (67 kDa) transmembrane protein (1). Immunolocalization studies indicate that **GLT1** is localized in astroglial cells throughout the brain and spinal cord.

**Source of Antigen and Antibodies**

<b>Antigen</b>	16-aa peptide from <b>rat GLT1 (1)</b> ; (Gene Accession#P31596) <b>Designation (# GLT11-P, control/blocking peptide)</b> conjugated to KLH. <b>Epitope location</b> ~ N-terminus Cytoplasmic domain
<b>Ab Host/type</b>	Rabbit, Polyclonal unpurified antiserum ( <b>#GLT11-S</b> ) and IgG, purified over antigen-agarose (Cat <b># GLT11-A</b> )
<b>2-Ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute** powder in 100 ul PBS

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder** in PBS at 1mg/ml

**Control/blocking peptide**

100 ug/100ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder** in PBS at 1 mg/ml.

**Storage**

**Short-term:** unopened, undiluted liquid vials at -200C and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** 1:1K-5K for antiserum and 1-10 ug/ml for affinity pure IgG using Chemiluminescence technique. This antibody has been referenced in many studies (2).

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies).

**Histochemistry & Immunofluorescence:** We recommend the use of affinity purified IgG at 2-10 ug/ml (see published refs 2 using this antibody).

**Specificity & Cross-reactivity**

The rat GLT11-P peptide is 100% conserved in mouse and 81% conserved in human. The GLT11-P has no significant sequence homology to other glutamate transporters. Antibody crossreactivity in various species is not established. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at the web site:

**General References: (1).** Pines, G. (1992) Nature 360:464-467; Kanner, BI (1993) FEBS Lett. 325, 95-99; Zhang, Y. (1994) J. Biol. Chem. 269:19573-19577; Grunewald, M. and Kanner, B. (1995) J. Biol. Chem. 270:17017-17024; Arriza, J. L. (1994) J. Neurosci 14:5559-5569; Krischner, M. A. (1994) Genomics 24:218-224;

**(2) Citations of ADI's Antibodies** (see web site for updated list)

Ueda Y, 2000, Exp. Brain Res. 133, 334-340, WB,,  
Voutsinos-Porche, 2003, Neuron, 37, 275-286, WB, IHC  
Rodriguez-Kern A, 2003, Neurochemistry Intl. 43, 363-370, WB,  
Yang Y, 2004, J. Neurochem.,; 88: 1361 - 1372., WB,,  
Ueda Y et al, 2001, J. Neurochem. 2001 76: 892-900, WB,,  
Levenson J, 2002, Nature Neuroscience 5, 155 - 161, WB,,  
Samuelsson C, 2000, Neurosci. lett. 289, 185-188, ,  
Li, Shuxin, 1999, J. Neurosci. 99, RC16, 1-9, , IHC,  
Ueda Y, 2002, Mol. brain. Res. 104, 120-126, WB,,  
Palos, Teresa P., 1999, J. Neurochemistry 73, 1012, WB,, IHC,

*\*This product is for In vitro research use only.*

**Related material available from ADI**

GLT11-S-A-P 71226S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)