

Product Specification Sheet

**Glucose Transporter 1 (Glut-1) Antibodies**

Cat. # GT13-P	Human Glut-1 Control/blocking Peptide	<b>SIZE:</b> 100 ug
Cat. # GT13-S	Rabbit Anti-Human Glut-1 antiserum	<b>SIZE:</b> 100 ul
Cat. # GT13-A	Rabbit Anti-Human Glut-1 IgG (affinity pure)	<b>SIZE:</b> 100 ug

Most mammalian cells transport glucose through a family of membrane proteins known as glucose transporters. Molecular cloning of these glucose transporters has identified a family of closely related genes that encodes at least 7 proteins (**Glut-1 to Glut-13**, Mol. Wt. 40-80 kDa) and Sodium glucose co-transporter-1 (SGLT-1, 662 amino acids; ~75 kDa). Individual member of this family have identical predicted secondary structures with 12 transmembrane domains. Both N and c-termini are predicted to be cytoplasmic. Most differences in sequence homology exist within the four hydrophilic domains that may play a role in tissue-specific targeting. Glut isoforms differ in their tissue expression, substrate specificity and kinetic characteristics. **Glut-1** mediates glucose transport into red cells, and throughout the blood brain barrier, and supply glucose to most cells.

**FUNCTION:** Facilitative glucose transporter. This isoform may be responsible for constitutive or basal glucose uptake. Has a very broad substrate specificity; can transport a wide range of aldoses including both pentoses and hexoses.

**SUBCELLULAR LOCATION:** Cell membrane; Multi-pass membrane protein (By similarity). Melanosome (By similarity). Note=Localizes primarily at the cell surface

**PTM:** Phosphorylated upon DNA damage, probably by ATM or ATR

**SIMILARITY:** Belongs to the major facilitator superfamily. Sugar transporter (TC 2.A.1.1) family. Glucose transporter subfamily Protein name Solute carrier family 2, facilitated glucose transporter member 1 ; Synonyms Glucose transporter type 1, erythrocyte/brain; GLUT-1, GT1, Glut1, Gene name : Slc2a1

**Source of Antigen and Antibodies**

<b>Antigen</b>	15-aa peptide <b>Human Glut-1 protein accession #P11166, refs 1; Designation (GT13-P, control peptide)</b> conjugated to KLH; <b>Epitope location</b> ~ C-terminal, Cytoplasmic domain
<b>Ab Host/type</b>	Rabbit, Polyclonal unpurified antiserum (#GT13-S) and IgG, purified over antigen-agarose (Cat # GT13-A)
<b>2-Ab</b>	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
<b>-ve control IgG</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**

100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute powder in 100 ul PBS**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer:** PBS+0.1% BSA  
**Reconstitute powder in PBS at 1mg/ml**

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at -200C and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Recommended Usage**

**Western Blotting** (1:1K-5K for antiserum and 1-10 ug/ml for affinity pure IgG). (see published papers refs 4)

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

**Histochemistry:** We recommend the use of affinity purified antibody at 2-10 ug/ml. (see published papers refs 4)

**Specificity & Cross-reactivity**

The human GT13 peptide is 100% conserved in pig, bovine, rabbit, mouse, and rat. Antibody crossreactivity in various species is not experimentally confirmed. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at: [www.4adi.com/data/abblock.html](http://www.4adi.com/data/abblock.html)).

**General References:** 1. Mueckler, M (1985) Science 229, 941-945; Haspel (1986) J. Biol. Chem. 263, 398-403; Birnbaum, (1986) 83, 5784-5788. . Piper (1991) Am. J. Physiol. 260, C570-C580; Harris . (1992) Proc. Natl. Acad. Sci. 89, 7556-7560.;

**2) Citations of ADI's antibodies for Glut-1**

Choeiri C, 2005, Neuroscience, 130, 591-600, WB,  
Majumder PK, 2004, Nature Medicine10, 594 - 601, WB,  
Mekhail K, 2004, Nature Cell Biology 6, 642 - 647, WB,  
Bindra RJ, 2004, Mol. Cell. Biol., Oct 2004; 24: 8504 - 8518, WB,,  
Lewis MD, 2004, Oncogene23, 2315 - 2323, WB,  
Gnudi L, 2003, Hypertension. 42(1):19-24, WB, IHC  
Gunaratnam L, 2003, J. Biol. Chem., 278: 44966 - 44974, WB,  
Greenberg CC, 2003, J. Biol. Chem.,278: 30835 - 30842., WB,,  
Hansen WJ eat al, 2002, Mol. Cell. Biol. 22: 1947-1960, WB,,  
Kondo K , 2002, Cancer Cell 1, 237-246, WB,  
Bindra RJ, 2002, Cancer Res. 62: 3014 - 3019, WB,,  
Tong, Haiyan, 2000, J. Biol. Chem. 2000 275: 11981-11986, w,  
Cockman, M E , 2000, J. Biol. Chem. 275: 25733-25741, WB,,  
Schoenfeld, Alan, 1998, PNAS 95: 8817-8822, WB,,

\*This product is for In vitro research use only.

Antibodies for Glut 1-13 & SGLT-1/2

GT13-S-A-P 70911A

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)