

Product Specification Sheet

**Glucose Transporter 7 (rGlut-7) Antibodies**

Cat. # GT71-P	Rat Glut-7 Control/blocking Peptide #1	<b>SIZE:</b> 100 ug
Cat. # GT71-A	Rabbit Anti-Rat Glut-7 IgG #1 (aff pure)	<b>SIZE:</b> 100 ug
Cat. # GT71-S	Rabbit Anti-Rat Glut-7 antiserum #1	<b>SIZE:</b> 100 ul

Most mammalian cells transport glucose through a family of membrane proteins known as glucose transporters. Molecular cloning of these glucose transporters has identified a family of closely related genes that encodes at least 7 proteins (**Glut-1 to Glut-13**, Mol. Wt. 40-80 kDa) and Sodium glucose co-transporter-1 (SGLT-1, 662 amino acids; ~75 kDa). Individual member of this family have identical predicted secondary structures with 12 transmembrane domains. Both N and C-termini are predicted to be cytoplasmic. Most differences in sequence homology exist within the four hydrophilic domains that may play a role in tissue-specific targeting.

Glut isoform differ in their tissue expression, substrate specificity and kinetic characteristics. **Glut-1** mediates glucose transport into red cells, and throughout the blood brain barrier, and supply glucose to most cells. **Glut-2** provides glucose to the liver and pancreatic cells. **Glut-3** is the main transporter in neurons, whereas **Glut-4** is primarily expressed in muscle and adipose tissue and regulated by insulin. **Glut-5** transports fructose in intestine and testis. **Glut-6/SLC2A6**, originally described as glut-9 has been reassigned as Glut-6. Glut-6 (human 507 aa; ~ 48 kDa; chromosome 9q34) is most closely related to Glut-8 (~45% homology). It is highly expressed in brain, spleen, and peripheral leukocytes. **Glut-7**, expressed in liver and other gluconeogenic tissues, mediates glucose flux across endoplasmic reticulum membrane. **Glut-8** is found in adult testis and placenta. Human **Glut-9** is expressed in spleen, peripheral leukocytes and brain. Human **Glut-10** (541 aa, chromosome 20q13.1; ~30-35% homology with Glut-3 and Glut-8) has been identified as a candidate gene for NIDDM susceptibility. It is widely expressed with highest levels in liver and pancreas. **Glut-11** (496 aa, chromosome 22q11.2; ~41% identity with Glut-5) is expressed in heart and skeletal muscle. **Glut-12** (human 617 aa, monkey 621 aa; ~ 50 kDa; ~30% homology with Glut-4 and 40% with Glut-10) is expressed in skeletal muscle, adipose tissue, and small intestine. **Glut-13** or proton myo-inositol transporter (**HMIT**; human 629 aa, rat 618 aa, ~75-90 kDa/67 kDa protein) is highly expressed in glial cells and some neurons. Glut-13 transport activity was specific for myo-inositol. Rat HMIT is ~35% identical to rat GlutX1.

**Source of Antigen and Antibodies**

<b>Antigen</b>	14-aa peptide from rat <b>Glut-7</b> ; (Gene Accession #S24344) <b>Designation (GT71-P, control peptide)</b> conjugated to KLH
<b>Location</b>	~C-terminus of rat Glut-7
<b>Ab Host/type</b>	Rabbit, Polyclonal; Unpurified antiserum (cat #GT71-S) Aff pure IgG (cat # <b>GT71-A</b> )
<b>2-ab</b>	<b>Goat Anti-rabbit IgG-HRP</b> cat # 20320 (AP, biotin, FITC conjugates also available)
<b>-ve control</b>	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

**Form & Storage of Antibodies/Peptide Control**

**Antiserum (unpurified)**  
100ul solution lyophilized powder  
Supplied in Buffer: 0.05% azide  
**Reconstitute powder in 100 ul PBS**

**Affinity pure IgG**

100 ug/100ul solution lyophilized powder  
Supplied in **Buffer: PBS+0.1% BSA**  
**Reconstitute powder in PBS at 1 mg/ml**

**Control/blocking peptide**

100 ug/100 ul solution lyophilized powder  
Supplied in Buffer: PBS pH 7.5,  
**Reconstitute powder in PBS at 1 mg/ml.**

**Storage**

**Short-term:** unopened, undiluted liquid vials at -20OC and powder at 4oC or -20oC..

**Long-term:** at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

**Stability:** 6-12 months at -20oC or below.

**Shipping:** 4oC for solutions and room temp for powder

**Recommended Usage**

**Western Blotting** (1:1K-5K for antiserum and 1-10 ug/ml for affinity pure IgG using Chemiluminescence technique).

**ELISA:** Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

**Histochemistry:** not tested.

**Specificity & Cross-reactivity**

Rat GT71 sequence is 64% conserved in human. Glut-7 and Glut-2 sequences are 68% identical in rats and show similar secondary and tertiary structures. However, Glut-7 contains an extra 6 amino acid (KKMKND) at the C-terminal end. This sequence is very similar to the last 6 amino acids of some of the microsomal UDP-glucuronosyltransferase isoenzymes (see refs. 1 and other refs. Cited therein). The KK-K- has been considered motif for the retention of transmembrane proteins (3) and is not found in Gluts 1-5. The crossreactivity of anti-rat GT71 antibody to other proteins with an endoplasmic reticulum motif and species crossreactivity has not been determined. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol at: [www.4adi.com/data/abblock.html](http://www.4adi.com/data/abblock.html)).

**General References:** 1. Waddell, I.D., et al (1992) Biochem. J. 286, 173-177 & (1991) Biochem. J. 275, 363-377; 3. Jackson, M.R., et al (1990) EMBO J., 9, 3153-3162

**Citations of ADI's antibodies for Glucose transporters** (see updated list at: [www.4adi.com/flr/glutsflr.html](http://www.4adi.com/flr/glutsflr.html))

\*This product is for In vitro research use only.

**Related material available from ADI**

Antibodies for Glut 1-13 & SGLT-1/2

GT71 S-A-P

Rev. 40128S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)