

Item Name: Human Adult Ovary Whole Tissue Lysate

Cat# HAL-1316

Size: 100 ug

Tissue type Ovary

Source Human

Buffer: Laemmli SDS-PAGE sample buffer (reduced containing 10mM DTT)

Store Total tissue protein lysate was prepared in Hepes buffer, pH 7.9, KCl, dexoycholate, NP-40%, EDTA, sucrose, and a cocktail of protease inhibitors. Protein was optimized and then it is mixed with 2X SDS-PAGE Laemmli sample buffer (reduced). Lysate is typically supplied at a concentration of 2 mg/ml (or see lot specific concn on the vial).. Store at 2-8°C for short term use. For extended storage, freeze working aliquots at -70°C. Repeated freezing and thawing is not recommended. Under proper storage conditions the shelf life is half a year from the date of receipt.

Recommended Usage

Western blot-Load 5-20 ug/lane for testing with various antibodies. Heat at 95oC for 2-5 minutes prior to loading on gels.

Lysate preparation

Tissue lysate was prepared by homogenization in modified RIPA buffer (150 mM sodium chloride, 50 mM Tris-HCl, pH 7.4, 1 mM ethylenediaminetetraacetic acid, 1 mM phenylmethylsulfonyl fluoride, 1% Triton X-100, 1% sodium deoxycholic acid, 0.1% sodium dodecylsulfate, 5 µg/ml of aprotinin, 5 µg/ml of leupeptin. Tissue and cell debris was removed by centrifugation. Protein concentration was determined with Bio-Rad protein assay. The product was boiled for 5 min in 1 x SDS sample buffer (50 mM Tris-HCl pH 6.8, 12.5% glycerol, 1% sodium dodecylsulfate, 0.01% bromophenol blue) containing 10mM DTT.

All items are for in vitro research use only.

Related items

| | |
|------------|---|
| #90100 | Western blot recycling kit |
| 80200 | Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate |
| 80201 | Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate |
| 80202 | Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate |
| 80215 | NuGlo ECL Substrate Kit for Western blot |
| 80400-100 | Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA) |
| 80400-100 | Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA) |
| #ACTB12-M | Mouse monoclonal anti-beta actin IgG |
| #G3PDH12-M | Mouse monoclonal anti-G3PDH (GAPDH) |

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HAL-1316 90713V

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com