

## Item Name Human Cell Line Array II

**Cat#** HCLA-17602


**Size:** 1 package of 5 slides

**Cell Lines** HL-60, Jurkat, K-562, MOLT-4, Raji and THP-1

**Source** Human

**Buffer:** In order to keep the antigens in their native forms, cells were fixed in 4% paraformaldehyde.

**Store** Can be stored for up to 3 months at 4°C under desiccated conditions.

<p><b>Human Cell Array II</b></p>  <p>HL-60 Jurkat K-562 MOLT-4 Raji THP-1</p>	<p><b>Recommended Usage</b></p> <p>This product is provided for immunocytochemical analysis of protein and protein-protein interaction.</p> <p><b>Cell Array preparation</b></p> <p>Human Cell Array II carries 6 different human cell lines, including HL-60, Jurkat, K-562, MOLT-4, Raji and THP-1. Cells are harvested at the log phase of growth. In order to keep the antigens in their native forms, cells were fixed in 4% paraformaldehyde. The cells were arrayed on a 12-well (5 mm) adhesive coated slide, with each well's surface specifically treated to enhance cellular attachment and to minimize background staining. Each well loaded with approximately <math>7.5 \times 10^3</math> cells. A follow-up heat dehydration process was done to ensure the attachment of cells and the stability of cellular proteins.</p> <p>All items are for in vitro research use only</p>
--	---

### Related items

#90100	Western blot recycling kit
80200	Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
80201	Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
80202	Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
80215	NuGlo ECL Substrate Kit for Western blot
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
#ACTB12-M	Mouse monoclonal anti-beta actin IgG
#G3PDH12-M	Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HCLA-17602 80703S

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)