

**Item Name** Raji Cell Slide

**Cat#** HCLS-17007


**Size:** 1 package of 5 slides

**Cell Lines** (11yrs, Male) B lymphocyte Burkitt's lymphoma

**Source** Human

**Buffer:** In order to keep the antigens in their native forms, cells were fixed in 4% paraformaldehyde.

**Store** Can be stored for up to 3 months at 4°C under desiccate conditions.

	<p><b>Recommended Usage</b></p> <p>This product is provided for immunocytochemical analysis of protein and protein-protein interaction.</p> <p><b>Cell Slide preparation</b></p> <p>Human Raji cells were cultured in RPMI 1640 medium with 2 mM L-glutamine and harvested at the log phase of growth. The Raji cells were fixed in 4% paraformaldehyde and arrayed on a 12-well (5 mm) adhesive coated slide, with each well's surface specifically treated to enhance cellular attachment and to minimize background staining. Each well was loaded with approximately <math>7.5 \times 10^3</math> Raji cells. A follow-up heat dehydration process was done to ensure the attachment of Raji cells and the stability of cellular proteins.</p> <p>All items are for in vitro research use only</p>
--	--

**Related items**

- #90100 Western blot recycling kit
- 80200 Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
- 80201 Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
- 80202 Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
- 80215 NuGlo ECL Substrate Kit for Western blot
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- #ACTB12-M Mouse monoclonal anti-beta actin IgG
- #G3PDH12-M Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HCLS-17007 80703S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)