

Product Specification Sheet

Alpha Defensin 1 (HDEFA1) Antibodies

Cat. # HDEFA11-P	HDEFA1 Alpha defensin 1 Control Peptide	SIZE: 100 ug
Cat. # HDEFA11-S	Rabbit Anti-HDEFA1 antiserum	SIZE: 100 ul
Cat. # HDEFA11-A	Rabbit Anti-HDEFA1 IgG (affinity pure)	SIZE: 100 ug

Defensins are a large family of broad-spectrum antimicrobial peptides, identified originally in leukocytes of rabbits and humans. Defensins, cationic/polar peptides (30-35 aa; 3-4 kDa), are distinguished by a conserved tri-disulfides and a largely β -Sheet structure. The genes encoding human α and β -defensins are clustered in a contiguous segment of chromosome 8p23. Defensins are classified into two families designated α - and β - based on distinctive, although similar, tri-disulfide linkages in the peptides. β -defensins are slightly larger and differ in the position and arrangement of 3 disulfides. In humans, six α -defensin (**cryptidins**), **HD 1-6** (HD1-4 are also known as **HNP1-4** for Human Neutrophil Peptides), and two β -defensins, **HBD-1** and **HBD-2**, have been identified.

α -defensins are encoded by DEFA1-6 genes, whereas human β -defensins are encoded by the DEFB1 and DEFB2 genes. HD1-4 are expressed in neutrophils, whereas HD5 and HD6 are expressed in epithelial cells of the intestinal and reproductive tract. HD1-3/HNP1-3 make up about 30% of the neutrophil's total granule protein. The HNPs are 29 to 30 amino acids long and are identical in sequence except at the N-terminal amino acid. **DEFA1** and **DEFA3** differ from each other only by their terminal ala and asp, respectively. **DEFA2** is identical to both except that it has 29 instead of 30 amino acids.

Source of Antigen and Antibodies

Antigen	Human alpha 1 defensins (HDEFA1) is synthesized from 94 precursor (mature peptide 65-94). A 16-aa peptide (designated HDEFA11-P; control peptides) Epitope location ~ C-terminus of human mature alpha 1 defensin-1 (1) conjugated to KLH
Ab Host/type	Rabbit, Polyclonal unpurified antiserum (#HDEFA11-S) and IgG, purified over antigen-agarose (Cat # HDEFA11-A)
2-Ab	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
-ve control	# 20009-1, Rabbit (non-immune) IgG, purified, suitable for ELISA, Western, IHC as -ve control

Form & Storage of Antibodies/Peptide Control

Antiserum (unpurified)

100ul solution lyophilized powder
Supplied in Buffer: 0.05% azide
Reconstitute powder in 100 ul PBS

Affinity pure IgG

100 ug/100ul solution lyophilized powder
Supplied in **Buffer:** PBS+0.1% BSA
Reconstitute powder in PBS at 1mg/ml

Control/blocking peptide

100 ug/100 ul solution lyophilized powder
Supplied in Buffer: PBS pH 7.5,
Reconstitute powder in PBS at 1 mg/ml.

Storage

Short-term: unopened, undiluted liquid vials at 20°C and powder at 4°C or -20°C..

Long-term: at -20°C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at -20°C or below.

Shipping: 4°C for solutions and room temp for powder

Recommended Usage

Western Blotting (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure antibody using ECL technique).

ELISA: Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

Histochemistry & Immunofluorescence: Not tested

Specificity & Cross-reactivity

The **HDEFA11-P** control peptides is common to human HNP1-3. The peptide has significant sequence homology with monkey alpha defensins 1-3. No significant homology exists with human alpha-defensins 4-6. The HDEFA11-P control peptide is available for antibody blocking to confirm specificity of antibodies. Antibody cross-reactivity in various species has not been studied.

General References: Selstead ME et al (1985) J Clin Invest. 76, 1436-1439; mars MW et al (1988) Blood 71, 1713-1719; Daher KA et al (1988) PNAS 85, 7327-7331; Linzmeier R et al (1993) FEBS lett. 321, 267-273; Wilson CL et al (1999) Science 286, 113-117; Mallow EB et al (1996) J. Biol. Chem. 271, 4038-4045; Ganz T and Lehrer RI (1998) Current Opinion Biol. 10, 41-44;.

*This product is for In vitro research use only.

Related material available from ADI

Antibodies alpha and beta-defensins and MMP7

Study distribution of proteins in pre-made **Kidney blots** from 7 defined regions of rat kidney

HDEFA11-S-A-P

71215S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com