

Item Name Human Tumor Tissue Blot IV

Cat# HTPB-1545

Size: 1 blot

Samples Human Tumor Tissue (kidney, liver, lymphoma, Non-Hodgkin's, spleen, thymoma, thyroid, tonsil)
Up to 8 custom samples can be run on a single blot

Protein Load: 15 ug total protein per lane

Gel type 4-20 % SDS-APGE Gradient gel (reducing) or Customized upon request

Mol wt markers 4-200 Kda (specified on the membrane)

Membrane type Nitrocellulose

Blocking **None;** Appropriate blocking agents (5% BSA or milk or other buffer should be used)

Store at 4oC in supplied bag

Recommended Usage

Human Tumor Tissue Blot IV is provided for Western Blot analysis. Before use in conventional immunoblotting assays, these blots need to be re-hydrated using any standard nitrocellulose membrane re-hydration procedure. For Western Blot, the blot can be blocked with 1% BSA or 5% nonfat dry milk for 1 hour. Incubation takes 30 minutes to two hours with primary antibody, and 30 minutes to one hour with secondary antibody at room temperature.

Blots may re-stripped with Western blot recycling kit # 901000 and probed with control antibodies (beta-actin, GAPDH) or other primary antibodies.

Notes: All human derived materials has been tested for HIV, hepatitis and found to be negative. Nevertheless all materials should be handled with precautions and considered hazardous. It should be used and disposed accordingly.

All items are for in vitro research use only.

Related items

#90100	Western blot recycling kit
80200	Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
80201	Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
80202	Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
80215	NuGlo ECL Substrate Kit for Western blot
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
#ACTB12-#ACTB12-M	Mouse monoclonal anti-beta actin IgG
#G3PDH12-M	Mouse monoclonal anti-G3PDH (GAPDH)

HTPB-1545 80702S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com