

Item Name Human Breast Slide (Poorly differentiated carcinoma)


Cat# HTS-10013

Size: 1 package of 5 slides

Tissue Lines Breast

Source Human

Store Stable for up to two years when stored in a cool, dry area.

	<p>Recommended Usage</p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p>Tissue Slide preparation</p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p>
------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Related items

#90100 Western blot recycling kit
 80200 Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
 80201 Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
 80202 Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
 80215 NuGlo ECL Substrate Kit for Western blot
 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
 #ACTB12-M Mouse monoclonal anti-beta actin IgG
 #G3PDH12-M Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-10013 80707S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
 Email: customerservice@lifetechindia.com Website: www.lifetechindia.com