

Item Name Human Colorectal Slide (normal)

Cat# HTS-10701

Size: 1 package of 5 slides

Tissue Lines Colorectal

Source Human

Store Stable for up to two years when stored in a cool, dry area.

| | |
|--|---|
|  | <p>Recommended Usage</p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p>Tissue Slide preparation</p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p> |
|--|---|

Related items

| | |
|------------|---|
| #90100 | Western blot recycling kit |
| 80200 | Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate |
| 80201 | Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate |
| 80202 | Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate |
| 80215 | NuGlo ECL Substrate Kit for Western blot |
| 80400-100 | Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA) |
| 80400-100 | Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA) |
| #ACTB12-M | Mouse monoclonal anti-beta actin IgG |
| #G3PDH12-M | Mouse monoclonal anti-G3PDH (GAPDH) |

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-10701 80707S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com