

Item Name Human Thyroid Slide (abnormal)


Cat# HTS-11018

Size: 1 package of 5 slides

Tissue Lines Thyroid

Source Human

Store Stable for up to two years when stored in a cool, dry area.

	<p>Recommended Usage</p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p>Tissue Slide preparation</p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p>
--	---

Related items

#90100	Western blot recycling kit
80200	Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
80201	Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
80202	Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
80215	NuGlo ECL Substrate Kit for Western blot
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
80400-100	Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
#ACTB12-M	Mouse monoclonal anti-beta actin IgG
#G3PDH12-M	Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-11018 80707S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com