

**Item Name** Human Salivary Gland Tissue Slide (Mucoepidermoid carcinoma)


**Cat#** HTS-12109

**Size:** 1 package of 5 slides

**Tissue Lines** Salivary Gland

**Source** Human

**Store** Stable for up to two years when stored in a cool, dry area.

	<p><b>Recommended Usage</b></p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p><b>Tissue Slide preparation</b></p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p>
------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

**Related items**

- #90100 Western blot recycling kit
- 80200 Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
- 80201 Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
- 80202 Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
- 80215 NuGlo ECL Substrate Kit for Western blot
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- #ACTB12-M Mouse monoclonal anti-beta actin IgG
- #G3PDH12-M Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-12109 80707S

**India Contact:**

**Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)