

Item Name Human Soft Tissue Slide (Benign)


Cat# HTS-12303

Size: 1 package of 5 slides

Tissue Lines Soft Tissue

Source Human

Store Stable for up to two years when stored in a cool, dry area.

	<p>Recommended Usage</p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p>Tissue Slide preparation</p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p>
------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Related items

- #90100 Western blot recycling kit
- 80200 Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate
- 80201 Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate
- 80202 Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate
- 80215 NuGlo ECL Substrate Kit for Western blot
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)
- #ACTB12-M Mouse monoclonal anti-beta actin IgG
- #G3PDH12-M Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-12303 80707S

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com