

## Item Name Human Skin Tissue Slide (Benign)

**Cat#** HTS-12702

**Size:** 1 package of 5 slides

**Tissue Lines** Skin

**Source** Human

**Store** Stable for up to two years when stored in a cool, dry area.

	<p><b>Recommended Usage</b></p> <p>Human tissue slides can be used with any antibody shown to be effective in formalin-fixed tissue using a standard immunostaining protocol.</p> <p>This product has not been evaluated for effectiveness in ISH or FISH applications</p> <p><b>Tissue Slide preparation</b></p> <p>The single tissue slides contain 4 µm sections of formalin-fixed, paraffin-embedded human tissue (1-2 cm x 1-2 cm) affixed to Snowcoat X-tra (Surgipath) slides. Tissue diagnosis was performed by certified pathologists.</p> <p>All items are for in vitro research use only</p>
------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

### Related items

#90100 Western blot recycling kit  
 80200 Western blot Kit for Rabbit Primary Antibodies, Chemilum. Substrate  
 80201 Western blot Kit for Goat Primary Antibodies, Chemilum. Substrate  
 80202 Western blot Kit for Mouse Primary Antibodies, Chemilum. Substrate  
 80215 NuGlo ECL Substrate Kit for Western blot  
 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)  
 80400-100 Bovine Serum Albumin (BSA, protease-free) for western (makes ~2-L buffer @5% BSA)  
 #ACTB12-M Mouse monoclonal anti-beta actin IgG  
 #G3PDH12-M Mouse monoclonal anti-G3PDH (GAPDH)

RIPA buffer, mol wt markers, Mouse, rat and human tissue lysates and pre-made protein blots are also available.

HTS-12702 80707S

### India Contact:

### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
 Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)