


## Y-27632

<b>Cat. No.:</b>	HY-10071
<b>CAS No.:</b>	146986-50-7
<b>Molecular Formula:</b>	C <sub>14</sub> H <sub>21</sub> N <sub>3</sub> O
<b>Molecular Weight:</b>	247.34
<b>Target:</b>	ROCK; Apoptosis
<b>Pathway:</b>	Cell Cycle/DNA Damage; Cytoskeleton; Stem Cell/Wnt; TGF-beta/Smad; Apoptosis
<b>Storage:</b>	-20°C, protect from light, stored under nitrogen * In solvent : -80°C, 6 months; -20°C, 1 month (protect from light, stored under nitrogen)



### SOLVENT & SOLUBILITY

#### In Vitro

DMSO : 50 mg/mL (202.15 mM; Need ultrasonic)  
H<sub>2</sub>O : 5 mg/mL (20.22 mM; ultrasonic and warming and heat to 60°C)

Concentration	Mass		
	1 mg	5 mg	10 mg
<b>1 mM</b>	4.0430 mL	20.2151 mL	40.4302 mL
<b>5 mM</b>	0.8086 mL	4.0430 mL	8.0860 mL
<b>10 mM</b>	0.4043 mL	2.0215 mL	4.0430 mL

Please refer to the solubility information to select the appropriate solvent.

#### In Vivo

- Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline  
Solubility: ≥ 2.5 mg/mL (10.11 mM); Clear solution
- Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline)  
Solubility: ≥ 2.5 mg/mL (10.11 mM); Clear solution
- Add each solvent one by one: 10% DMSO >> 90% corn oil  
Solubility: ≥ 2.5 mg/mL (10.11 mM); Clear solution
- Add each solvent one by one: Saline  
Solubility: 2 mg/mL (8.09 mM); Clear solution; Need ultrasonic and warming and heat to 60°C

### BIOLOGICAL ACTIVITY

#### Description

Y-27632 is an orally active, ATP-competitive inhibitor of ROCK-I and ROCK-II, with K<sub>i</sub>s of 220 and 300 nM, respectively. Y-27632 attenuates Doxorubicin-induced apoptosis of human cardiac stem cells. Y-27632 also suppresses dissociation-induced apoptosis of murine prostate stem/progenitor cells. Y-27632 primes human induced pluripotent stem cells (hiPSCs) to selectively differentiate towards mesendodermal lineage via epithelial-mesenchymal transition-like modulation<sup>[1][2][3][4][5][6][7]</sup>.

<b>IC<sub>50</sub> &amp; Target</b>	ROCK-I 220 nM (Ki)	ROCK-II 300 nM (Ki)	PKN 3.1 μM (Ki)	Citron kinase 5.3 μM (Ki)
	PKCα 73 μM (Ki)	PKA 25 μM (Ki)		
<b>In Vitro</b>	<p>Y-27632 inhibits the ROCK family of kinases 100 times more potently than other kinases including protein kinase C, cAMP-dependent kinase and myosin light chain kinase. Y-27632 prolongs the lag time and delays the appearance of BrdU-labeled cells in a concentration-dependent manner, delays of about 1 and 4 h are noticed in the Swiss 3T3 cells treated with 10 and 100 μM Y-27632, respectively<sup>[1]</sup>.</p> <p>Y-27632 promotes neuronal differentiation of adipose tissue-derived stem cells (ADSCs). Compared to 1.0 and 2.5 μM Y-27632 induced groups, percentages of neuroal-like cells achieved a peak in the 5.0 μM Y-27632 induced group<sup>[2]</sup>.</p> <p>Extracellular matrix (ECM) molecules decreases apoptosis markers and inhibiting the ROCK pathway blocks ECM stimulated actin cortical mat reformation and increases apoptosis in embryonic corneal epithelial cells<sup>[8]</sup>.</p> <p>MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>			
<b>In Vivo</b>	<p>Y-27632 (5 and 10 mg/kg) significantly prolongs the onset time of myoclonic jerks when compare with saline group. Y-27632 (5 and 10 mg/kg) significantly prolongs the onset time of clonic convulsions when compare with saline group<sup>[3]</sup>. Treatment with Dimethylnitrosamine (DMN) causes a significant decrease in rat body and liver weight (DMN-S group) compared with control animals (S-S group). Oral Y27632 (30 mg/kg) essentially prevents this DMN-induced rat body and liver weight loss (DMN-Y group)<sup>[4]</sup>.</p> <p>MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>			

## PROTOCOL

### Kinase Assay <sup>[1]</sup>

Recombinant ROCK-I, ROCK-II, PKN, or citron kinase is expressed in HeLa cells as Myc-tagged proteins by transfection using Lipofectamine, and is precipitated from the cell lysates by the use of 9E10 monoclonal anti-Myc antibody coupled to G protein-Sepharose. Recovered immunocomplexes are incubated with various concentrations of [<sup>32</sup>P]ATP and 10 mg of histone type 2 as substrates in the absence or presence of various concentrations of either Y-27632 or Y-30141 at 30°C for 30 min in a total volume of 30 μL of the kinase buffer containing 50 mM HEPES-NaOH, pH 7.4, 10 mM MgCl<sub>2</sub>, 5 mM MnCl<sub>2</sub>, 0.02% Brij 35, and 2 mM dithiothreitol. PKCα is incubated with 5 μM [<sup>32</sup>P]ATP and 200 μg/mL histone type 2 as substrates in the absence or presence of various concentrations of either Y-27632 or Y-30141 at 30°C for 10 min in a kinase buffer containing 50 mM Tris-HCl, pH 7.5, 0.5 mM CaCl<sub>2</sub>, 5 mM magnesium acetate, 25 μg/mL phosphatidyl serine, 50 ng/mL 12-O-tetradecanoylphorbol-13-acetate and 0.001% leupeptin in a total volume of 30 μL. Incubation is terminated by the addition of 10 μL of 43 Laemmli sample buffer. After boiling for 5 min, the mixture is subjected to SDS-polyacrylamide gel electrophoresis on a 16% gel. The gel is stained with Coomassie Brilliant Blue, and then dried. The bands corresponding to histone type 2 are excised, and the radioactivity is measured<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### Cell Assay <sup>[1]</sup>

HeLa cells are plated at a density of 3×10<sup>4</sup> cells per 3.5-cm dish. The cells are cultured in DMEM containing 10% FBS in the presence of 10 mM Thymidine for 16 h. After the cells are washed with DMEM containing 10% FBS, they are cultured for an additional 8 h, and then 40 ng/mL of Nocodazole is added. After 11.5 h of the Nocodazole treatment, various concentrations of Y-27632 (0-300 μM), Y-30141, or vehicle is added and the cells are incubated for another 30 min<sup>[1]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

### Animal Administration <sup>[3][4]</sup>

Mice<sup>[3]</sup>

Male, inbred Swiss albino mice (2-3 months old) weighing 25-30 g are used. Mice are injected with a sub-convulsive dose of PTZ (35 mg/kg, i.p.) (on Mondays, Wednesdays and Fridays) of each week for a total of 11 injections. After each PTZ injection, mice are observed for 30 min and the occurrence of convulsive activity is recorded. After 30 min, the mice are then injected with either Fasudil (25 mg/kg, i.p.) or Y-27632 (5 mg/kg, i.p.) and returned to their home cages until the next injection. Control mice for Fasudil and Y-27632 receives saline.

Rats<sup>[4]</sup>

Male Wistar Kind A rats (200-250 g) are used. DMN (1 g/mL) is diluted ten times with saline (final concentration 1%) and 10 mg/kg per day of DMN is injected intraperitoneally (i.p.) on the first 3 days of each week for 4 weeks. Y27632 is given orally once per day at a dose of 30 mg/kg for 4 weeks starting on the day of the first injection of DMN. The dose of 30 mg/kg corrects hypertension in several rat models without toxicity. Twenty rats are randomized into four experimental groups (n=5 in each group) as follows: (1) S-S (injection of saline i.p. and oral administration of saline); (2) S-Y (injection of saline i.p. and oral administration of Y27632); (3) DMN-S (DMN i.p. and oral administration of saline); (4) DMN-Y (DMN i.p. and oral administration of Y27632). The rats are weighed every week. They are sacrificed at the end of the fourth week and the liver is excised. In addition, a blood sample is taken immediately before the rats are sacrificed.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## CUSTOMER VALIDATION

- Nature. 2021 Dec 2.
- Science. 2020 Dec 4;370(6521):eaay2002.
- Cell Stem Cell. 2020 Oct 1;27(4):646-662.e7.
- Cell Stem Cell. 2020 Jul 2;27(1):125-136.e7.
- Cell Stem Cell. 2020 Jun 4;26(6):845-861.e12.

See more customer validations on [www.MedChemExpress.com](http://www.MedChemExpress.com)

## REFERENCES

- [1]. Ishizaki T, et al. Pharmacological properties of Y-27632, a specific inhibitor of rho-associated kinases. *Mol Pharmacol.* 2000 May;57(5):976-83.
- [2]. Xue ZW, et al. Rho-associated coiled kinase inhibitor Y-27632 promotes neuronal-like differentiation of adult human adipose tissue-derived stem cells. *Chin Med J (Engl).* 2012 Sep;125(18):3332-5.
- [3]. Inan S, et al. Antiepileptic effects of two Rho-kinase inhibitors, Y-27632 and fasudil, in mice. *Br J Pharmacol.* 2008 Sep;155(1):44-51.
- [4]. Tada S, et al. A selective ROCK inhibitor, Y27632, prevents dimethylnitrosamine-induced hepatic fibrosis in rats. *J Hepatol.* 2001 Apr;34(4):529-36.
- [5]. Maldonado M, et al. ROCK inhibitor primes human induced pluripotent stem cells to selectively differentiate towards mesendodermal lineage via epithelial-mesenchymal transition-like modulation. *Stem Cell Res.* 2016 Sep;17(2):222-227.
- [6]. Kan L, et al. Rho-Associated Kinase Inhibitor (Y-27632) Attenuates Doxorubicin-Induced Apoptosis of Human Cardiac Stem Cells. *PLoS One.* 2015;10(12):e0144513. Published 2015 Dec 8.
- [7]. Zhang L, et al. ROCK inhibitor Y-27632 suppresses dissociation-induced apoptosis of murine prostate stem/progenitor cells and increases their cloning efficiency. *PLoS One.* 2011;6(3):e18271. Published 2011 Mar 28.
- [8]. Svoboda KK, et al. ROCK inhibitor (Y27632) increases apoptosis and disrupts the actin cortical mat in embryonic avian corneal epithelium. *Dev Dyn.* 2004;229(3):579-590.

**Caution: Product has not been fully validated for medical applications. For research use only.**

**India Contact:**

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)