


## LY3009120

<b>Cat. No.:</b>	HY-12558		
<b>CAS No.:</b>	1454682-72-4		
<b>Molecular Formula:</b>	C <sub>23</sub> H <sub>29</sub> FN <sub>6</sub> O		
<b>Molecular Weight:</b>	424.51		
<b>Target:</b>	Raf; Autophagy		
<b>Pathway:</b>	MAPK/ERK Pathway; Autophagy		
<b>Storage:</b>	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month



### SOLVENT & SOLUBILITY

#### In Vitro

DMSO : ≥ 38 mg/mL (89.51 mM)  
 \* "≥" means soluble, but saturation unknown.

Preparing Stock Solutions	Solvent Concentration	Mass		
		1 mg	5 mg	10 mg
	1 mM	2.3557 mL	11.7783 mL	23.5566 mL
	5 mM	0.4711 mL	2.3557 mL	4.7113 mL
	10 mM	0.2356 mL	1.1778 mL	2.3557 mL

Please refer to the solubility information to select the appropriate solvent.

#### In Vivo

1. Add each solvent one by one: 5% DMSO >> 95% (20% SBE-β-CD in saline)  
 Solubility: 2.5 mg/mL (5.89 mM); Suspended solution; Need ultrasonic

### BIOLOGICAL ACTIVITY

#### Description

LY3009120 is a pan RAF inhibitor which inhibits BRAF<sup>V600E</sup>, BRAF<sup>WT</sup> and CRAF<sup>WT</sup> with IC<sub>50</sub>s of 5.8, 9.1 and 15 nM, respectively.

#### IC<sub>50</sub> & Target

BRAF <sup>V600E</sup> 5.8 nM (IC <sub>50</sub> )	Braf 9.1 nM (IC <sub>50</sub> )	CRAF 15 nM (IC <sub>50</sub> )
---	------------------------------------	-----------------------------------

#### In Vitro

In the whole-cell based KiNativ assay, LY3009120 shows affinity to each RAF isoform with the IC<sub>50</sub> of 44, 31-47 and 42 nM for ARAF, BRAF and CRAF respectively. LY3009120 exhibits anti-proliferative effects on cell lines harboring BRAF<sup>V600E</sup>, KRAS<sup>G13</sup> and KRAS<sup>G12</sup> mutations. LY3009120 (1 μM) inhibits the phosphorylation of both MEK1/2 and ERK1/2 in cell lines with high basal levels of pMEK1/2 and pERK1/2 (RKO and HCT 116)<sup>[1]</sup>. LY3009120 shows inhibitory effect on tumor cells such as BxPC-3, NCI-H2405 and OV-90 cell lines. LY3009120 (0.01 μM) demonstrates potent and dose-dependent inhibition of phospho-MEK and ERK in all three cell lines. LY3009120 demonstrates a concentration-dependent cell growth inhibition with IC<sub>50</sub>

values of 0.04, 0.087, and 0.007  $\mu\text{M}$  against H2405, BxPC-3, and OV-90 cells, respectively<sup>[2]</sup>. LY3009120 inhibits BRAF<sup>WT</sup>, CRAF<sup>WT</sup>, BRAF<sup>V600E</sup>, and BRAF<sup>V600E+G468A</sup> with the IC<sub>50</sub> values of 9.1, 15, 5.8, and 17 nM, respectively. LY3009120 induces BRAF-CRAF dimerization but inhibits the phosphorylation of downstream MEK and ERK. LY3009120 also inhibits various forms of RAF dimers including BRAF or CRAF homodimers<sup>[3]</sup>. LY3009120 gives only very minor activation at very low doses, with near complete inhibition of phospho-ERK at concentrations above 100 nM<sup>[4]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

#### In Vivo

LY3009120 (20 mg/kg bid) displays significant activity in in vivo BRAF<sup>mut</sup> and KRAS<sup>mut</sup> CRC xenograft models. In Colo 205 xenografts (BRAF<sup>mut</sup>), LY3009120 results in statistically significant tumor regression, while treatment of HCT 116 xenografts (KRAS<sup>mut</sup>) results in statistically significant inhibition of tumor growth. LY3009120 treatment reduces pMEK1/2 in all HT-29 xenografts and reduces pERK1/2 in the majority of HT-29 xenografts<sup>[1]</sup>. LY3009120 (15 or 30 mg/kg) achieves almost complete tumor growth regression, and inhibits downstream phospho-MEK and ERK by approximately 70% and 60%, respectively, in the H2405 model<sup>[2]</sup>.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## PROTOCOL

#### Cell Assay <sup>[4]</sup>

Briefly, cells are grown in McCoy's 5A supplemented with 10% characterized fetal bovine serum at 37°C, 5% CO<sub>2</sub>, and 95% humidity. Cells are allowed to expand until 75-90% confluency at which point they are subcultured or harvested for assay use. A serial dilution of test compound is dispensed into a 384-well black clear bottom plate in triplicate. Six-hundred-twenty-five cells are added per well in 50  $\mu\text{L}$  of complete growth medium in the 384-well plate. Plates are incubated for 67 h at 37°C, 5% CO<sub>2</sub>, and 95% humidity. At the end of the incubation period, 10  $\mu\text{L}$  of a 440  $\mu\text{M}$  solution of resazurin in PBS is added to each well of the plate and plates are incubated for an additional 5 h at 37°C, 5% CO<sub>2</sub>, and 95% humidity. Plates are read on a Synergy2 reader using an excitation of 540 nm and an emission of 600 nm. Data are analyzed using Prism software to calculate IC<sub>50</sub> values.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

#### Animal Administration <sup>[2]</sup>

Briefly,  $5 \times 10^6$  to  $10 \times 10^6$  tumor cells in a 1:1 Matrigel mix (0.2 mL total volume) are injected subcutaneously into the right hind flank of female NIH nude rats. After tumors reach a desired size of approximately 300 mm<sup>3</sup>, animals are randomized into groups of 8 for efficacy studies. Drugs (LY3009120 or PLX4032) are administered orally (gavage) in 0.6-mL volume of vehicle with the dose schedules. Tumor growth and body weight are monitored over time to evaluate efficacy and signs of toxicity.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## CUSTOMER VALIDATION

- Cell Syst. 2020 Nov 18;11(5):478-494.e9.
- Oncogene. 2018 Oct;37(43):5719-5734.
- Clin Sci (Lond). 2019 Apr 16;133(8):919-932.
- Cancer Sci. 2018 Jan;109(1):121-131.
- bioRxiv. 2019 Sep.

See more customer validations on [www.MedChemExpress.com](http://www.MedChemExpress.com)

## REFERENCES

[1]. Vakana E, et al. LY3009120, a panRAF inhibitor, has significant anti-tumor activity in BRAF and KRAS mutant preclinical models of colorectal cancer. Oncotarget. 2017

[2]. Chen SH, et al. Oncogenic BRAF Deletions That Function as Homodimers and Are Sensitive to Inhibition by RAF Dimer Inhibitor LY3009120. *Cancer Discov.* 2016 Mar;6(3):300-15

[3]. Peng SB, et al. Inhibition of RAF Isoforms and Active Dimers by LY3009120 Leads to Anti-tumor Activities in RAS or BRAF Mutant Cancers. *Cancer Cell.* 2015 Sep 14;28(3):384-98

[4]. Henry JR, et al. Discovery of 1-(3,3-dimethylbutyl)-3-(2-fluoro-4-methyl-5-(7-methyl-2-(methylamino)pyrido[2,3-d]pyrimidin-6-yl)phenyl)urea (LY3009120) as a pan-RAF inhibitor with minimal paradoxical activation and activity against BRAF or RAS mutant tumors

---

**Caution: Product has not been fully validated for medical applications. For research use only.**

**India Contact:**

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)