


Amsacrine hydrochloride

Cat. No.:	HY-13551A
CAS No.:	54301-15-4
Molecular Formula:	C ₂₁ H ₂₀ ClN ₃ O ₃ S
Molecular Weight:	429.92
Target:	Topoisomerase; Autophagy
Pathway:	Cell Cycle/DNA Damage; Autophagy
Storage:	4°C, protect from light
	* In solvent : -80°C, 6 months; -20°C, 1 month (protect from light)



H-Cl

SOLVENT & SOLUBILITY

In Vitro	DMSO : 62.5 mg/mL (145.38 mM; Need ultrasonic)			
		Solvent Concentration	Mass	
			1 mg	5 mg
			10 mg	
Preparing Stock Solutions	1 mM	2.3260 mL	11.6301 mL	23.2601 mL
	5 mM	0.4652 mL	2.3260 mL	4.6520 mL
	10 mM	0.2326 mL	1.1630 mL	2.3260 mL
Please refer to the solubility information to select the appropriate solvent.				
In Vivo	1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.08 mg/mL (4.84 mM); Clear solution			

BIOLOGICAL ACTIVITY

Description	Amsacrine hydrochloride (m-AMSA hydrochloride; acridinyl anisidine hydrochloride) is an inhibitor of topoisomerase II , and acts as an antineoplastic agent which can intercalates into the DNA of tumor cells.
IC ₅₀ & Target	Topoisomerase II
In Vitro	Amsacrine (mAMSA) blocks HERG currents in HEK 293 cells and Xenopus oocytes in a concentration-dependent manner, with IC ₅₀ values of 209.4 nM and 2.0 μM, respectively. Amsacrine (mAMSA) causes a negative shift in the voltage dependence of both activation (-7.6 mV) and inactivation (-7.6 mV). HERG current block by Amsacrine (mAMSA) is not frequency dependent ^[1] . In vitro studies of normal human lymphocytes with various concentrations of Amsacrine (mAMSA), show both increased levels of chromosomal aberrations, ranging from 8% to 100%, and increase SCEs, ranging from 1.5 times the normal at the lowest concentration studied (0.005 μg/mL) to 12 times the normal (0.25 μg/mL) ^[3] . Amsacrine (mAMSA)-induced apoptosis of U937 cells is characterized by caspase-9 and caspase-3 activation, increased intracellular Ca ²⁺ concentration, mitochondrial depolarization, and MCL1 down-

	<p>regulation. Amsacrine induces MCL1 down-regulation by decreasing its stability. Further, amsacrine-treated U937 cells show AKT degradation and Ca²⁺-mediated ERK inactivation^[4].</p> <p>MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>
In Vivo	<p>In animals treated with different doses of amsacrine (0.5-12 mg/kg), the frequencies of micronucleated polychromatic erythrocytes increase significantly after treatment with 9 and 12 mg/kg. Furthermore, the present study demonstrates for the first time that Amsacrine (mAMSA) has high incidences of clastogenicity and low incidences of aneugenicity whereas nocodazole has high incidences of aneugenicity and low incidences of clastogenicity during mitotic phases in vivo^[2].</p> <p>MCE has not independently confirmed the accuracy of these methods. They are for reference only.</p>

PROTOCOL

Animal Administration ^[2]

Mice^[2]

Amsacrine (mAMSA) is investigated in three separated experiments. In the first experiment, animals are treated by intraperitoneal injection with 0.5, 1.5 and 4.5 mg/kg of Amsacrine (mAMSA) and bone marrow is sampled 24 h after treatment. Preliminary negative MN results at this sampling time lead to the use of 30 h sampling time for Amsacrine (mAMSA). Thus, in the second experiment, mice are treated with 0.5, 1.5 and 4.5 mg/kg of Amsacrine (mAMSA) and bone marrow is sampled 30 h after treatment. The doses and sampling times for Amsacrine are chosen by reference to earlier studies and the selected doses are within the dose range used for human chemotherapy. The results again show that the micronuclei frequency in the bone marrow of mice is not affected by treatment with any of the selected doses of the test agent, at 30 h sampling time, thus, in the third experiment, mice are treated with 6, 9 and 12 mg/kg of Amsacrine and bone marrow is sampled 24 and 30 h after treatment.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- J Mol Med (Berl). 2019 Aug;97(8):1183-1193.

See more customer validations on www.MedChemExpress.com

REFERENCES

- [1]. Thomas D, et al. Inhibition of cardiac HERG currents by the DNA topoisomerase II inhibitor amsacrine: mode of action. Br J Pharmacol. 2004 Jun;142(3):485-94.
- [2]. Attia SM. Molecular cytogenetic evaluation of the mechanism of genotoxic potential of amsacrine and nocodazole in mouse bone marrow cells. J Appl Toxicol. 2013 Jun;33(6):426-33.
- [3]. Kao-Shan CS, et al. Cytogenetic effects of amsacrine on human lymphocytes in vivo and in vitro. Cancer Treat Rep. 1984 Jul-Aug;68(7-8):989-97.
- [4]. Lee YC, et al. Amsacrine-induced apoptosis of human leukemia U937 cells is mediated by the inhibition of AKT- and ERK-induced stabilization of MCL1. Apoptosis. 2017 Mar;22(3):406-420.

Caution: Product has not been fully validated for medical applications. For research use only.

Tel: 609-228-6898

Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA