


SN-38

Cat. No.:	HY-13704		
CAS No.:	86639-52-3		
Molecular Formula:	C ₂₂ H ₂₀ N ₂ O ₅		
Molecular Weight:	392.4		
Target:	Topoisomerase; ADC Cytotoxin; Autophagy		
Pathway:	Cell Cycle/DNA Damage; Antibody-drug Conjugate/ADC Related; Autophagy		
Storage:	Powder	-20°C	3 years
		4°C	2 years
	In solvent	-80°C	6 months
		-20°C	1 month



SOLVENT & SOLUBILITY

In Vitro

DMSO : 25 mg/mL (63.71 mM; Need ultrasonic)
 H₂O : < 0.1 mg/mL (insoluble)

Preparing Stock Solutions	Solvent Concentration	Mass		
		1 mg	5 mg	10 mg
	1 mM	2.5484 mL	12.7421 mL	25.4842 mL
	5 mM	0.5097 mL	2.5484 mL	5.0968 mL
	10 mM	0.2548 mL	1.2742 mL	2.5484 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline
 Solubility: 2.5 mg/mL (6.37 mM); Suspended solution; Need ultrasonic
- Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline)
 Solubility: 2.5 mg/mL (6.37 mM); Suspended solution; Need ultrasonic
- Add each solvent one by one: 10% DMSO >> 90% corn oil
 Solubility: 2.5 mg/mL (6.37 mM); Suspended solution; Need ultrasonic

BIOLOGICAL ACTIVITY

Description

SN-38 (NK012) is an active metabolite of the Topoisomerase I inhibitor Irinotecan. SN-38 (NK012) inhibits DNA and RNA synthesis with IC₅₀s of 0.077 and 1.3 μM, respectively.

IC₅₀ & Target

Topoisomerase I Camptothecins

In Vitro

The IC₅₀ values for LoVo, HCT116, and HT29 cell lines is 20 nM, 50 nM, 130 nM, respectively. In all three SN-38 (NK012)

resistant cell lines Top1 activity is maintained in the presence of high concentrations of SN-38^[2].
MCE has not independently confirmed the accuracy of these methods. They are for reference only.

In Vivo

SN-38 (NK012), the active and toxic metabolite of the anticancer prodrug Irinotecan. At 30 minutes after administration, Irinotecan plasma concentrations in Slco1a1b(-/-) mice are 1.9-fold higher than in the wild-type mice (1.89 vs. 1.01 μ M, respectively), whereas SN-38 (NK012) plasma concentrations of Slco1a1b(-/-) mice are 8-fold higher compare with wild-type mice (0.4 μ g/mL vs. 0.05 μ g/mL, respectively). Overall plasma exposure [AUC₍₅₋₂₄₀₎] of Irinotecan is 1.7-fold higher in Oatp1a1b knockout mice versus wild-type mice (209.8 \pm 6.7 vs. 120.9 \pm 4.4 μ M/min; P<0.01), and 2.9-fold higher for SN-38 (50 \pm 2.9 vs. 12 \pm 2 μ M/min; P<0.001)^[3].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Cell Assay ^[2]

In vitro SN-38 (NK012) sensitivity is determined using the MTT assay. Cells are seeded in 96-well plates, and a range of SN-38 (NK012) concentrations is added the following day. Following 48 h of drug exposure, the medium is discarded and the plates are incubated with medium containing MTT (0.5 mg/mL) for 3 h. Acidified (0.02 M HCl) sodium dodecyl sulphate (20 %) is added to dissolve the formed formazan. Optical density at 570 nm (and 670 nm for background) is measured, and the cell viability is calculated in percent compared to untreated cells. Experiments are repeated three times and the mean IC₅₀ value \pm standard deviation is determined. Relative resistance for each resistant cell line is calculated by dividing the mean IC₅₀ value of the resistant cell line by the mean IC₅₀ value of the corresponding parental cell line^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration ^[3]

Mice^[3]

Female wild-type, Slco1a1b(-/-)(Oatp1a1b knockout), Slco1a1b(-/-);1B1(tg), and Slco1a1b(-/-);1B3(tg) (liver-specific OATP1B1 and OATP1B3 humanized transgenic) mice of comparable genetic background (>99% FVB) between 8 and 14 weeks of age are used. Irinotecan (20 mg/mL in water-based solution containing NaOH, lactic acid, and sorbitol) is diluted with saline (to 2 mg/mL) for administration of 10 mg/kg; 5 μ L/g bodyweight are administered intravenously to mice. SN-38 (NK012) is dissolved in DMSO (1 mg/mL) and 1 μ L/g body weight is administered intravenously to mice to achieve a dosage of 1 mg/kg. The experiments are terminated by isoflurane anaesthesia, heparin-blood sampling by cardiac puncture followed by cervical dislocation and tissue collection. Blood samples are centrifuged at 5,200 \times g for 5 minutes at 4°C and plasma is collected and stored at -30°C until analysis.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- J Control Release. 2020 Jul 21;326:387-395.
- Carbohydr Polym. 2020 May 1;235:115983.
- Int J Nanomedicine. 2020, 15: 6839-6854.
- INT J PHARMACEUT. 2020 May.
- Int J Pharm. 2019 Oct 5;569:118588.

See more customer validations on www.MedChemExpress.com

REFERENCES

[1]. Wallin A, et al. Anticancer effect of SN-38 on colon cancer cell lines with different metastatic potential. Oncol Rep. 2008 Jun;19(6):1493-8.

[2]. Jensen NF, et al. Characterization of DNA topoisomerase I in three SN-38 resistant human colon cancer cell lines reveals a newpair of resistance-associated mutations.

Caution: Product has not been fully validated for medical applications. For research use only.

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com