


Ethambutol dihydrochloride

Cat. No.:	HY-B0535A
CAS No.:	1070-11-7
Molecular Formula:	C ₁₀ H ₂₆ Cl ₂ N ₂ O ₂
Molecular Weight:	277.23
Target:	Bacterial; Antibiotic
Pathway:	Anti-infection
Storage:	4°C, sealed storage, away from moisture * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture)



SOLVENT & SOLUBILITY

In Vitro

DMSO : 100 mg/mL (360.71 mM; Need ultrasonic)
H₂O : ≥ 50 mg/mL (180.36 mM)
* "≥" means soluble, but saturation unknown.

Preparing Stock Solutions	Solvent Concentration	Mass	1 mg	5 mg	10 mg
		1 mM	3.6071 mL	18.0356 mL	36.0711 mL
	5 mM	0.7214 mL	3.6071 mL	7.2142 mL	
	10 mM	0.3607 mL	1.8036 mL	3.6071 mL	

Please refer to the solubility information to select the appropriate solvent.

In Vivo

1. Add each solvent one by one: PBS
Solubility: 100 mg/mL (360.71 mM); Clear solution; Need ultrasonic

BIOLOGICAL ACTIVITY

Description

Ethambutol dihydrochloride (Emb dihydrochloride) is a bacteriostatic antimycobacterial agent, which obstructs the formation of cell wall by inhibiting arabinosyl transferases. Target: Antibacterial Ethambutol dihydrochloride (Emb dihydrochloride) directly affects two polymers, arabinogalactan (AG) and lipoarabinomannan (LAM) in *Mycobacterium smegmatis*. In *M. smegmatis*, Ethambutol inhibits synthesis of arabinan completely and inhibits AG synthesis most likely as a consequence of this; more than 50% of the cell arabinan is released from the bacteria following Ethambutol treatment, whereas no galactan is released. Ethambutol main targets against embB gene product in *M. avium*. Ethambutol induces 60% changes in the embB gene in *M. tuberculosis* resistant mutants [1]. Ethambutol dihydrochloride (Emb dihydrochloride) is effective against actively growing microorganisms of the genus *Mycobacterium*, including *M. tuberculosis*. Nearly all strains of *M. tuberculosis* and *M. kansasii* as well as a number of strains of the *M. avium* complex (MAC) are sensitive to Ethambutol. [1] Ethambutol dihydrochloride (Emb dihydrochloride) is potency against *M. tuberculosis* (H37Rv) with MIC of 0.5 µg/mL in vitro [2]. Ethambutol is efficient on treatment of mycobacterial-infected macrophages. When *M. tuberculosis* infected macrophages are treated with 6 µg/mL Ethambutol, the log CFUs following treatment for 3 days is 4.17, while value in

control group is 4.8. The MICs for *M. avium* (MTCC 1723) and *M. smegmatis* (MTCC 6) are 15 µg/mL and 0.18 µg/mL, respectively. Ethambutol is efficient in animal model. 100 mg/kg Ethambutol given orally 15 days post i.v. infection 1 ×/week for 5 weeks, induces a lower log CFU compared with untreated (4.59 vs 5.07) [3].

CUSTOMER VALIDATION

- ACS Chem Biol. 2021 Dec 15.

See more customer validations on www.MedChemExpress.com

REFERENCES

- [1]. Ethambutol. *Tuberculosis (Edinb)*, 2008. 88(2): p. 102-5.
- [2]. Rastogi, N., V. Labrousse, and K.S. Goh, In vitro activities of fourteen antimicrobial agents against drug susceptible and resistant clinical isolates of *Mycobacterium tuberculosis* and comparative intracellular activities against the virulent H37Rv strain in human macrophages. *Curr Microbiol*, 1996. 33(3): p. 167-75.
- [3]. Kaur, D. and G.K. Khuller, In vitro, ex-vivo and in vivo activities of ethambutol and sparfloxacin alone and in combination against mycobacteria. *Int J Antimicrob Agents*, 2001. 17(1): p. 51-5.
-

Caution: Product has not been fully validated for medical applications. For research use only.

India Contact:
Life Technologies (India) Pvt. Ltd.
306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com