

## β-Amyloid (1-42), human TFA

<b>Cat. No.:</b>	HY-P1363	
<b>CAS No.:</b>	107761-42-2	
<b>Molecular Formula:</b>	C <sub>205</sub> H <sub>312</sub> F <sub>3</sub> N <sub>55</sub> O <sub>62</sub> S	
<b>Molecular Weight:</b>	4628.06	
<b>Sequence:</b>	Asp-Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala-Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly-Val-Val-Ile-Ala	
<b>Sequence Shortening:</b>	DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVVIA	
<b>Target:</b>	Amyloid-β	
<b>Pathway:</b>	Neuronal Signaling	
<b>Storage:</b>	Powder    -80°C    2 years -20°C    1 year In solvent   -80°C    6 months -20°C    1 month	

DAEFRHDSGYEVHHQKLVFFAEDVGSNKGAIIGLMVGGVVIA (TFA salt)

### SOLVENT & SOLUBILITY

#### In Vitro

DMSO : 50 mg/mL (10.80 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Concentration	Mass		
		1 mg	5 mg	10 mg
	1 mM	0.2161 mL	1.0804 mL	2.1607 mL
	5 mM	0.0432 mL	0.2161 mL	0.4321 mL
	10 mM	0.0216 mL	0.1080 mL	0.2161 mL

Please refer to the solubility information to select the appropriate solvent.

### BIOLOGICAL ACTIVITY

#### Description

β-Amyloid (1-42), human TFA (Amyloid β-Peptide (1-42) (human) TFA) is a 42-amino acid peptide which plays a key role in the pathogenesis of Alzheimer disease<sup>[1]</sup>.

#### In Vitro

β-Amyloid Aggregation Guidelines (Following is our recommended protocol. This protocol only provides a guideline, and should be modified according to your specific needs).

1. Solid Aβ peptide was dissolved in cold hexafluoro-2-propanol (HFIP). The peptide was incubated at room temperature for at least 1h to establish monomerization and randomization of structure.
2. The HFIP was removed by evaporation, and the resulting peptide was stored as a film at -20 or -80 °C.
3. The resulting film was dissolved in anhydrous DMSO at 5 mM and then diluted into the appropriate concentration and buffer (serum- and phenol red-free culture medium) with vortexing.
4. Next, the solution was age 48h at 4-8 °C. The sample was then centrifuged at 14000g for 10 min at 4-8 °C; the soluble oligomers were in the supernatant. The supernatant was diluted 10-200-fold for experiments.

---

Methods vary depends on the downstream applications.  
MCE has not independently confirmed the accuracy of these methods. They are for reference only.

## CUSTOMER VALIDATION

- Yonsei Med J. 2019 Jul;60(7):640-650.

See more customer validations on [www.MedChemExpress.com](http://www.MedChemExpress.com)

## REFERENCES

- [1]. Solntseva EI, et al. Impact of amyloid- $\beta$  peptide (1-42) on voltage-gated ion currents in molluscan neurons. Bull Exp Biol Med. 2011 Oct;151(6):671-4.
  - [2]. Barucker C, et al. Nuclear translocation uncovers the amyloid peptide A $\beta$ 42 as a regulator of gene transcription. J Biol Chem. 2014 Jul 18;289(29):20182-91.
  - [3]. Stefania Sabella, et al. Capillary electrophoresis studies on the aggregation process of beta-amyloid 1-42 and 1-40 peptides. Electrophoresis. 2004 Oct;25(18-19):3186-94.
- 

**Caution: Product has not been fully validated for medical applications. For research use only.**

India Contact:  
Life Technologies (India) Pvt. Ltd.  
306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)