

ImmunoStrip® for Potato virus X (PVX)



Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

The PVX ImmunoStrip® detects all strains of PVX in potato leaf tissue.

The test is highly sensitive (comparable to ELISA methods) and can be completed in as little as a few minutes with no equipment or expertise.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB4 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB4 extraction buffer purchased separately

Option

ImmunoStrip® for Potato virus X (PVX) 5 strips & buffer filled bags
ISK 10000/0005

ImmunoStrip® for Potato virus X (PVX), 25 strips & buffer filled bags
ISK 10000/0025

ImmunoStrips® only for Potato virus X (PVX) ***BUFFER FILLED BAGS NOT INCLUDED***
STX 10000/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com