

ImmunoStrip® for Pepino mosaic virus (PepMV)



The *Pepino mosaic virus* (PepMV) ImmunoStrip® is used to detect the presence of PepMV in solanaceous crops. PepMV is a member of the Potexvirus genus known for their non-enveloped, flexuous, filamentous-shaped virus particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

	<p><i>Included in ISK Kit:</i></p> <ul style="list-style-type: none"> • ImmunoStrips® • Sample bags containing SEB1 buffer • User Guide <p><i>Included in STX Item:</i></p> <ul style="list-style-type: none"> • ImmunoStrips® • User Guide <p><input type="checkbox"/> Option</p>
<input type="checkbox"/>	ImmunoStrip® for Pepino mosaic virus (PepMV), 5 strips & buffer filled bags ISK 13001/0005
<input type="checkbox"/>	ImmunoStrip® for Pepino mosaic virus (PepMV), 25 strips & buffer filled bags ISK 13001/0025
<input type="checkbox"/>	ImmunoStrips® only for Pepino mosaic virus (PepMV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 13001/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com