

## Agdia ImmunoStrip® for Potato virus Y (PVY)

The *Potato virus Y* (PVY) ImmunoStrip® is used to detect the presence of PVY in potato leaves and sprouts. PVY is a member of the Potyvirus genus known for their non-enveloped, flexuous, filamentous virus particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

If you would prefer to add a field-deployable molecular assay to your testing toolkit, please check out Agdia's [PVY AmplifyRP® XRT](#) method.

### *Included in ISK Kit:*

- ImmunoStrips®
- Sample bags containing [SEB4 buffer](#)
- [User Guide](#)

### *Included in STX Item:*

- ImmunoStrips®
- [User Guide](#)

Option
ImmunoStrip® for Potato virus Y (PVY), 5 strips and buffer filled bags ISK 20001/0005
ImmunoStrip® for Potato virus Y (PVY), 25 strips and buffer filled bags ISK 20001/0025
ImmunoStrips® only for Potato virus Y (PVY) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 20001/0025

### India Contact:

#### **Life Technologies (India) Pvt. Ltd.**

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)