

ImmunoStrip® for Potato virus X & Potato virus Y (PVX & PVY)

The *Potato virus X* and *Potato virus Y* (PVX & PVY) dual pathogen ImmunoStrip® is used to detect the presence of PVX and PVY in potato leaves. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

If you would prefer to add a field-deployable molecular assay to your testing toolkit, please check out Agdia's [PVY AmplifyRP® XRT](#) method.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB4 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for Potato virus X & Potato virus Y (PVX & PVY), 5 strips and buffer filled bags ISK 41300/0005
<input type="checkbox"/>	ImmunoStrip® for Potato virus X & Potato virus Y (PVX & PVY), 25 strips and buffer filled bags ISK 41300/0025
<input type="checkbox"/>	ImmunoStrips® only for Potato virus X & Potato virus Y (PVX & PVY) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 41300/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com