

Agdia ImmunoStrip® for Lettuce mosaic virus (LMV)

The *Lettuce mosaic virus* (LMV) ImmunoStrip® is used to detect the presence of LMV in lettuce leaf and seed. LMV is a member of the Potyvirus genus known for their non-enveloped, flexuous, filamentous virus particles. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

<input type="checkbox"/>	Option
<input type="checkbox"/>	ImmunoStrip® for LMV, 5 strips and buffer filled bags ISK 45801/0005
<input type="checkbox"/>	ImmunoStrip® for LMV, 25 strips and buffer filled bags ISK 45801/0025
<input type="checkbox"/>	ImmunoStrip® for LMV ***BUFFER FILLED BAGS NOT INCLUDED*** STX 45801/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com