

## Agdia ImmunoStrip® for *Pantoea stewartii* (Pstew)

The *Pantoea stewartii* (Pstew) ImmunoStrip® is used to detect the presence of Pstew, the causal agent of Stewart's Wilt, in corn, jackfruit, and wheat leaves. This test is not validated to test corn seed. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

### Option

ImmunoStrip® for Pstew, 5 strips and buffer filled bags  
ISK 52000/0005

ImmunoStrip® for Pstew, 25 strips and buffer filled bags  
ISK 52000/0025

ImmunoStrip® for Pstew \*\*\*BUFFER FILLED BAGS NOT INCLUDED\*\*\*  
STX 52000/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)