

## ImmunoStrip® for Iris yellow spot virus (IYSV)

The *Iris yellow spot virus* (IYSV) ImmunoStrip® is used to detect the presence of IYSV in allium species. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

### Options

Option
ImmunoStrip® for Iris yellow spot virus (IYSV), 5 strips & buffer filled bags ISK 60500/0005
ImmunoStrip® for Iris yellow spot virus (IYSV), 25 strips & buffer filled bags ISK 60500/0025
ImmunoStrips® only for Iris yellow spot virus (IYSV) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 60500/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)