

ImmunoStrip® for Tobacco ringspot virus (TRSV)



Test Label: NA
Test Format: ImmunoStrip®
Capture Reagent: Polyclonal
Detection Reagent: Polyclonal

Intended Use:

The TRSV ImmunoStrip® is an on-site tool to quickly identify plants infected with *Tobacco ringspot virus* (TRSV).

TRSV is a member of the genus *Nepovirus* and can infect a wide range of ornamental, fruit, vegetable and field crops. According to the European and Mediterranean Plant Protection Organization (EPPO) ([Bulletin](#)) TRSV can be transmitted mechanically and vectored by at least one nematode species and a variety of insects. In addition it can be transmitted by seed in several plant species.

The TRSV ImmunoStrip® has been tested against more than 75 susceptible hosts with no observable non-specific reactions. When tested against infected tissue sources this assay demonstrated equivalent diagnostic sensitivity and similar analytical sensitivity to that of our commercially available TRSV ELISA.

The TRSV ImmunoStrip® test is an ideal tool for on-site screening of plants. Use it in scouting regimes or to test plants before purchasing or accepting a shipment of plants.

Kit (ISK) Includes:

- ImmunoStrips®
- Sample bags containing SEB1 buffer
- [User Guide](#)

ImmunoStrips® Only (STX) Includes:

- ImmunoStrips®
- SEB1 extraction buffer purchased separately

Option

ImmunoStrip for TRSV, 5 strips & buffer filled bags
ISK 64001/0005

ImmunoStrip for TRSV, 25 strips & buffer filled bags
ISK 64001/0025

ImmunoStrips® only for Tobacco ringspot virus (TRSV) ***BUFFER FILLED
BAGS NOT INCLUDED***
STX 64001/0025

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444
Email: customerservice@lifetechindia.com Website: www.lifetechindia.com