

## ImmunoStrip® for *Xanthomonas axonopodis* pv. *citri* (Xac)

The *Xanthomonas axonopodis* pv. *citri* (Xac) ImmunoStrip® is used to detect the presence of Xac in citrus fruit. The Xac ImmunoStrips® can detect the Asiatic strain (strain A) of Xac. It does not detect strain A\*, strain A<sup>w</sup> (Wellington), or type-A etrog. ImmunoStrips® are the perfect screening tool for use in the field, greenhouse, and the lab.

### Included in ISK Kit:

- ImmunoStrips®
- Sample bags containing [SEB1 buffer](#)
- [User Guide](#)

### Included in STX Item:

- ImmunoStrips®
- [User Guide](#)

Option
ImmunoStrip® for <i>Xanthomonas axonopodis</i> pv. <i>citri</i> (Xac), 5 strips & buffer filled bags ISK 92200/0005
ImmunoStrip® for <i>Xanthomonas axonopodis</i> pv. <i>citri</i> (Xac), 25 strips & buffer filled bags ISK 92200/0025
ImmunoStrips® only for <i>Xanthomonas axonopodis</i> pv. <i>citri</i> (Xac) ***BUFFER FILLED BAGS NOT INCLUDED*** STX 92200/0025

### India Contact:

#### Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (India). Ph: +91-11-42208000, 42208111, 42208222 Mobile: 91-11-42208444  
Email: [customerservice@lifetechindia.com](mailto:customerservice@lifetechindia.com) Website: [www.lifetechindia.com](http://www.lifetechindia.com)