

Product Specification Sheet

KST1 Antibodies

Cat # KST11-P	Human KST1 control peptide # 1	SIZE: 100 ug
Cat # KST11-S	Rabbit Anti-Human KST1 Antiserum # 1	SIZE: 100 ul
Cat # KST11-A	Rabbit Anti-Human KST1 IgG # 1, aff pure	SIZE: 100 ug

The kidneys play a major role in the regulation of glucose levels. Kidneys filter approx. 180 g of glucose per day from the blood, and this is mostly reabsorbed back into the blood in the proximal tubules. Typically, glucose is first absorbed within epithelium by a specific transporter protein, Sodium glucose co transporters (SGLT), in the brush-border membrane and then it is transported out of the cell across the basolateral membranes by a facilitated sugar transporter (GLUTs). At least 3 members of SGLTs (SGLT1-3) have been cloned and characterized from various species. Individual member of this family have identical predicted secondary structures with up to 14 transmembrane domains. SGLT1-3 genes code for protein of approx 659-672 residues (calculated size of ~75 kDa). Both N and C-termini are predicted to be extracellular. There is approx 60-70% homology between SGLT1-3. SGLTs transport α-methyl-D-glucoside (α-MDG), a non-metabolized model substrate, in Na-dependent manner. SGLT1 does not discriminate α-MDG, glucose, and galactose. SGLT2/3 do not transport D-galactose efficiently.

A human gene, **KST1**, encoding a new member of SGLT family has been identified. KST1 (rabbit 674 aa, human 675 aa chromosome 16p12-p11), an ortholog of rabbit kidney SGLT or rKST, is expressed in brain, heart, muscle, kidney, liver and placenta. Due to its wide tissue expression and chromosomal location, mutations in KST1 have been implicated in ICCA (infantile convulsion and choreoathetosis) and BFIC (benign familial infantile convulsions).

Source of Antigen and Antibodies

Antigen	A 15-aa peptide sequence (designated KST11-P; control peptide) from the N-terminal, cytoplasmic domain of Human KST1 (1) was synthesized, coupled to KLH
Ab Host/type	Rabbit, Polyclonal antiserum # KST11-S and IgG, purified over antigen-agarose (Cat # KST11-A)
2-Ab	Cat # 20320, goat anti-rabbit IgG-HRP (AP, biotin, FITC conjugates also available).
-ve control IgG	Cat # 20009-1, Rabbit (non-immune) Serum IgG, purified, suitable for ELISA, Western, IHC as -ve control

Form & Storage of Antibodies/Peptide Control

Antiserum (unpurified)
100ul solution lyophilized powder
Supplied 0.05% azide, **Reconstitute** powder in 100 ul PBS

Affinity pure IgG

100 ug/100ul solution lyophilized powder
Supplied in **Buffer:** PBS+0.1% BSA
Reconstitute powder in PBS at 1mg/ml

Control/blocking peptide

100 ug/100 ul solution lyophilized powder
Supplied in Buffer: PBS pH 7.5,
Reconstitute powder in PBS at 1 mg/ml.

Storage

Short-term: unopened, undiluted liquid vials at -20OC and powder at 4oC or -20oC..

Long-term: at -20C or below in suitable aliquots after reconstitution. Do not freeze and thaw and store working, diluted solutions.

Stability: 6-12 months at -20oC or below.

Shipping: 4oC for solutions and room temp for powder

Recommended Usage

Western Blotting (1:1K-5K for neat serum and 1-10 ug/ml for affinity pure antibody using Chemiluminescence technique).

ELISA: Control peptide can be used to coat ELISA plates at 1 ug/ml and detected with antibodies (1:10-50K for neat serum and 0.5-1 ug/ml for affinity pure).

Histochemistry & Immunofluorescence: not tested.

Specificity & Cross-reactivity

The Human KST11-P peptide sequence is 66% conserved rabbit KST1. It has not yet been cloned from other species. No significant sequence homology exists with other SGLTs. Control peptide, because of its low mol. Wt (<3 kDa), is not suitable for Western. It should be used for ELISA or antibody blocking experiments (use 5-10 ug control peptide per 1 ug of aff pure IgG or 1 ul antiserum) to confirm antibody specificity (see detailed protocol see detailed protocol at the web site).

General References: (1) Roll P et al (2002) Gene 285, 141-148; Wright E (2001) Am. J. Physiol. Renal Physiol. 280, F10-F18 (review)

*This product is for In vitro research use only.

Antibodies to SGLT1-3, RS11, and Glut1-13

KST11-S-A-P 71213A

India Contact:

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi – 110034 (INDIA). Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com Website: www.lifetechindia.com