

Mouse Monoclonal Antibody (IgG) against SARS-CoV-2 Nucleocapsid protein

Catalog Number: LT362A14

Size: 1 mg

Type: mouse IgG1

Production:

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with recombinant SARS-CoV-2 Nucleoprotein

Formulation:

10 mM PB + 150 mM NaCl + 0.1% Proclin300, pH7.4.

Application/Usage:

Immunoassays

Preliminary pair recommendations in sandwich immunoassays (capture antibody)

Purification:

Protein A chromatography

Storage:

Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Material safety note:

This product is sold for research or further manufacturing use only. Standard Laboratory Practices should be followed when handling this material.

Life Technologies (India) Pvt. Ltd.

306, Aggarwal City Mall, Opposite M2K Pitampura, Delhi - 110034 (INDIA).

Ph: +91-11-42208000, 42208111, 42208222, Mobile: +91-9810521400, Fax: +91-11-42208444

Email: customerservice@lifetechindia.com

Website: www.lifetechindia.com