

## **Anti-CD16 Molecule (CD16) antibody (FITC)**

**Cat No:** LT3768742NIABU

**Antigen:** CD16 Molecule (CD16)

**Synonyms:** Fcgr3, Fcgr3a, CD16, FCGR3, FCGR3@, CD16b, FCG3, FCR-10, FCRIII, FCRIIIb, FcRIII

**Reactivity:** Cow (Bovine), Horse (Equine), Human, Pig (Porcine), Sheep (Ovine)

**Host: Mouse Clonality:** Monoclonal

**Conjugate: FITC Application:** Flow Cytometry (FACS)

**Product Details Clone: KD1 Isotype:** IgG2a **Cross-Reactivity:** Cow (Bovine), Horse (Equine), Pig (Porcine), Sheep (Ovine)

**Characteristics:** Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC)

**Target Details Alternative Name:** CD16

**Synonyms: FcRIII Gene ID:** 2214, 2215, 281766

**Application Notes:** Optimal working dilution should be determined by the investigator.

**Restrictions:** For Research Use only

**Handling Format:** Liquid

**Concentration:** 0.1 mg/ml