



NEST Polyclonal Antibody

Cat #: LT54495PBAKBA Size: 100µl

Product Information

	Product Name: NEST Polyclonal Antibody		
	Applications: WB, ELISA		Isotype: Rabbit IgG
	Reactivity: Human, Mouse, Rat		
REF	Catalog Number: LT54495PBAKBA	LOT	Lot Number: Refer to product label
	Formulation: Liquid		Concentration: 1 mg/ml
	Storage: Store at -20°C. Avoid repeated freeze / thaw cycles.		Note: Contain sodium azide.

Background: NES encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells. NES (Nestin) is a Protein Coding gene. Diseases associated with NES include Brain Stem Cancer and Ependymoblastoma. Among its related pathways are Embryonic and Induced Pluripotent Stem Cell Differentiation Pathways and Lineage-specific Markers and Neural Stem Cell Differentiation Pathways and Lineage-specific Markers. GO annotations related to this gene include structural molecule activity and intermediate filament binding. An important paralog of this gene is SYNM.

Application Notes: Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB (1:500-1:2000), ELISA (1:5000-1:20000).

Storage Buffer: PBS, pH 7.4, containing 0.02% Sodium Azide as preservative and 50% Glycerol.

Storage Instructions: Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.

Note: The product listed herein is for research use only and is not intended for use in human or clinical diagnosis. Suggested applications of our products are not recommendations to use our products in violation of any patent or as a license. We cannot be responsible for patent infringements or other violations that may occur with the use of this product.